



# lexCap CM 6FF lexCap DEAE 6FF

# lexCap SP 6FF lexCap Q 6FF

## Index

1. Product Description.....	1
2. Purification Procedure.....	2
3. Cleaning-in-Place.....	3
4. Troubleshooting.....	3
5. Related Products.....	3

## 1. Product Description

**CM, SP, DEAE, Q Beads 6FF** are part of Ion exchange Resin which is widely used in biomedical and bioengineering for separation and purification of proteins, nucleic acids and polypeptides. The base matrix of **CM, SP, DEAE, Q Beads 6FF** is 6% highly cross-linked agarose which gives the ion exchange resin chemical and physical stability. The characteristics such as capacity, elution behavior and pressure/flow rate are unaffected by the solutions commonly used in process chromatography and cleaning procedures, for details see table under each respective ion exchange resin. Fig.1 is the pressure/flow rate curves of the matrix.




Fig.1. A typical pressure/flow rate curve for **CM, SP, DEAE, Q Beads 6FF**

**IEXCap 6FF** is a prepacked ready to use column. It is packed with 1ml and 5ml of IEX medium. **IEXCap 6FF** can be adapted to all kinds of chromatography system, such as ÄKTA. It is easy to operate.

### CM Beads 6FF

**CM Beads 6FF** is a weak cation exchange resin. The ion exchange group is a carboxy methyl group, see below.




Table 1. Characteristics of **CM Beads 6FF**

Item	Description
Matrix	Highly cross-linked 6% agarose
Ion exchange type	Weak cation
Total ionic capacity	0.09-0.13mmol H <sup>+</sup> /ml medium
Particle Size	45-165 μm
Flow rate	300-600 cm/h
pH stability	4-13
Storage buffer	1×PBS containing 20% ethanol
Storage	4°C - 30°C

### SP Beads 6FF

**SP Beads 6FF** is a strong cation exchange resin. The ion exchange group is a sulphopropyl group, see below.






Table 2. Characteristics of SP Beads 6FF

Item	Description
Matrix	Highly cross-linked 6% agarose
Ion exchange type	Strong cation
Ion exchange capacity	0.18-0.25mmol H <sup>+</sup> /ml medium
Particle Size	45-165 μm
Flow rate	400-700 cm/h
pH stability	4-13
Storage buffer	20% ethanol, 0.2M sodium acetate
Storage	4°C - 30°C

### DEAE Beads 6FF

DEAE Beads 6FF is a weak anion exchange resin. The ion exchange group is a diethylaminoethyl group, see below.




Table 3. Characteristics of DEAE Beads 6FF

Item	Description
Matrix	Highly cross-linked 6% agarose
Ion exchange type	Weak anion
Ion exchange capacity	0.11-0.16mmol Cl <sup>-</sup> /ml medium
Particle Size	45-165 μm
Flow rate	300-600 cm/h
pH stability	2-12
Storage buffer	1× PBS containing 20% ethanol
Storage	4°C - 30°C

### Q Beads 6FF

Q Beads 6FF is a strong anion exchange resin. The ion exchange group is a quaternary amine group, see below.




Table 4. Characteristics of Q Beads 6FF

Item	Description
Matrix	Highly cross-linked 6% agarose
Ion exchange type	Strong anion
Ion exchange capacity	0.18-0.25mmol Cl <sup>-</sup> /ml medium
Particle Size	45-165 μm
Flow rate	400-700 cm/h
pH stability	2-12
Storage buffer	1×PBS containing 20% ethanol
Storage	4°C - 30°C

## 2. Purification Procedure

### 2.1 Buffer Preparation

Water and chemicals used for buffer preparation should be high purity. It is recommended to filter the buffers by passing them through a 0.22μm or 0.45μm filter before use.

### 2.2 Sample Preparation

It is recommended to filter the sample solution by passing them through a 0.22μm or 0.45μm filter before use.

### 2.3 Sample Purification

IEX 6FF is a prepacked, ready to use column. The prepacked column provides fast, simple and easy separations in a convenient format.

1) Fill the pump tubing with binding buffer. Connect the column to purification system, “drop to drop” to avoid introducing air into the column.





2) Wash the column with 10 column volumes of binding buffer.

3) Apply the sample, using a syringe fitted to the connector or by pumping it onto the column.

4) Wash with 5 to 10 column volumes of binding buffer or until no material appears in the effluent.

5) Elute with Elution Buffer using a stepwise or linear gradient. For one-step elution, 5 column volumes are usually enough. Other volumes may be required if the interaction is difficult to break. Linear gradient elution can be used to separate proteins of different binding strengths with a small gradient, such as 20 column volumes or more.

#### 2.4 Analysis

Identify the fractions containing the target protein. using UV absorbance, SDS-PAGE, or western blot.

### 3. Clean-in-Place

After each separation, elute any reversibly bound material either with a high ionic strength solution (e.g. 1M NaCl in buffer ) or by increasing pH. Regenerate the beads by washing with at least 5 bed volumes of buffer, or until the column effluent shows stable conductivity and pH values.

Cleaning-in-place (CIP) is a cleaning procedure that removes contaminants such as lipids, precipitates, or denatured proteins that may remain in the packed column after regeneration. Regular CIP also prevents the build-up of these contaminants in the media bed and helps to maintain the capacity, flow properties and general performance of the media.

A specific CIP protocol should be designed for each process according to the type of contaminants present. CIP cycle is generally recommended every 1-5 separation cycles.

#### Remove the ionically bound proteins

Wash with 3-4 column volumes of 2M NaCl. Contact time 10-15min.

#### Remove the precipitation or hydrophobically bound proteins or lipoproteins

Wash with at least 2 column volumes of 1M NaOH . Contact time 1-2h.

#### Remove lipids and very hydrophobic proteins

Wash with 2-4 column volumes of 0.5% non-ionic detergent, 70% ethanol or 30% isopropanol. Contact time 1-2h.

### 4. Troubleshooting

Problem	Probable Cause	Solution
Back pressure is too high	Column is clogged	Cleaning in place(part 3).
	Sample solution contains precipitate	Filtering the sample solution by passing them through a 0.22µm or 0.45µm filter.
Eluate is not pure	The medium repeat too much times.	Cleaning in place(part 3).
	Wash is not enough.	Increase the volume of Wash Buffer.

### 5. Related Products

Product	Cat. No.	Size
lexCap Q 6FF	SI001C11	1X1 ml
	SI001C51	5X1 ml
	SI001C15	1X5 ml
	SI001C55	5X5 ml
lexCap SP 6FF	SI003C11	1X1 ml
	SI003C51	5X1 ml
	SI003C15	1X5 ml
	SI003C55	5X5 ml
lexCap DEAE 6FF	SI005C11	1X1 ml
	SI005C51	5X1 ml
	SI005C15	1X5 ml
	SI005C55	5X5 ml
lexCap CM 6FF	SI007C11	1X1 ml
	SI007C51	5X1 ml





	SI007C15	1X5 ml
	SI007C55	5X5 ml

