

# PRODUCT BROCHURE

**Professional Manufacturer**  
of Chromatography Resins



# COMPANY PROFILE

Changzhou Smart-Lifesciences Biotechnology Co., Ltd., established in 2013 and affiliated with Changzhou Boyi Bio Group, is a high-tech enterprise focusing on the development and production of high-quality biological purification products. Our main offerings include chromatography resins, automated purification instruments, and supporting consumables. Adhering to the "Quality First" principle, the company has passed ISO9001 quality management system certification. Its mature production processes and standardized quality inspection systems ensure data reliability and stable quality, providing a guarantee for the safety and stability of product applications, and meeting the batch and quality requirements of industrial customers. The company's advanced technology also enables us to provide customers with professional technical support and process development services.



**3,000+**(L)

Maximum Single Batch



**34+**(PEC)

Patents



**2,500+**(CNT)

Partner Biopharmaceutical Enterprises



**20+**(YRS)

Development Experience

**Our products have been widely used in the following fields:**



Antibody Drugs



Polypeptide



Recombinant Protein



Gene Therapy



Blood Products



Vaccine



API



## STABLE SUPPLY

Single-batch capacity up to ton-scale ensures continuous and reliable supply chain operations.



## QUALITY ASSURANCE

Certified under ISO9001:2015 quality management system, ensuring process safety and regulatory compliance.



## EXPERT TEAM

Capable of meeting diverse microsphere customization needs, providing technical support for process development and validation.



## PROCESS SUPPORT

Dedicated after-sales service team committed to safeguarding your production processes.

**PROFESSIONAL  
MANUFACTURER OF  
CHROMATOGRAPHY RESINS**

# QUALITY SYSTEM CERTIFICATE



# CONTENTS |

## 01 Polysaccharide Chromatography Resin

### Affinity Chromatography Resin

|   |    |
|---|----|
| Histidine-tagged protein affinity purification resin      | 02 |
| Other metal ion chelating affinity purification resin     | 05 |
| GST fusion protein affinity purification resin            | 07 |
| MBP fusion protein affinity purification resin            | 08 |
| Biotin-tagged protein affinity purification resin         | 09 |
| Strep-tag II fusion protein affinity purification resin   | 10 |
| Flag, HA-c-Myc-tagged protein affinity purification resin | 12 |
| GFP, RFP-, YFP-tagged protein affinity purification resin | 14 |
| Protein A protein affinity purification resin             | 16 |
| Protein G, A/G, L protein affinity purification resin     | 19 |
| Pre-activated resin                                       | 21 |
| Other affinity purification resin                         | 23 |

### Ion Exchange Chromatography Resin

|   |    |
|---|----|
| Strong cation exchange chromatography resin | 28 |
| Weak cation exchange chromatography resin   | 29 |
| Strong anion exchange chromatography resin  | 29 |
| Weak anion exchange chromatography resin    | 30 |

### Hydrophobic Interaction Chromatography Resin

|  |    |
|--|----|
| Butyl, Octyl, Phenyl and other ligands for hydrophobic interaction chromatography (HIC) resins | 30 |
|--|----|

### Size Exclusion Chromatography Resin

|                                     |    |
|-------------------------------------|----|
| Agarose gel filtration resin        | 33 |
| Dextran-based gel filtration resins | 34 |
| Agarose-dextran composite SEC resin | 35 |
| Pre-packed desalting columns        | 35 |

### Multimodal Chromatography purification resin

|  |    |
|--|----|
|  | 37 |
|--|----|

## 02 Polymer-Based Chromatography Resin

|                           |    |
|---------------------------|----|
| Phmac series HIC resin    | 41 |
| Phmac series IEC resin    | 43 |
| Smart PS series IEC resin | 45 |
| Smart PS series RPC resin | 46 |

## 03 Other Chromatography Resin

|                      |    |
|----------------------|----|
| Smart C8/C18         | 49 |
| Monosilica C8/C18    | 50 |
| Hydroxyapatite (CHT) | 50 |


## 04 Magnetic Beads

|   |    |
|---|----|
| Magarose beads for protein purification | 53 |
| Magpoly beads for immunoprecipitation   | 54 |


## 05 Instruments and Consumables

|  |    |
|--|----|
| High-throughput purification equipment | 56 |
| Protein chromatography system          | 56 |
| Gravity Empty Column                   | 57 |
| Single-layer glass column              | 58 |
| Double-layer column (SXX series)       | 58 |

## Chromatography Resin Selection Guide



## Chromatography Resin Selection Guide



# Polysaccharide Chromatography Resin

Affinity Chromatography Resin

Ion Exchange Chromatography Resin

Hydrophobic Interaction Chromatography Resin

Size Exclusion Chromatography Resin

Multimodal Chromatography purification resin

# 01

## 1.1 Affinity Chromatography Resin

### Histidine-tagged protein affinity purification resin

Immobilized metal ion affinity chromatography (IMAC) is widely used in the purification of proteins and peptides. This technique relies on the specific interaction between transition metal ions and certain amino acid chains, such as histidine, tryptophan, and L-cysteine, that are exposed on the protein surface. IMAC media are particularly effective in separating proteins containing histidine residues. Key features of IMAC include high binding capacity, excellent stability, mild separation conditions, and the ability to quickly and reliably scale up the purification of His-tagged proteins. In recent years, IMAC has gained popularity for purifying histidine-tagged proteins expressed in various systems, including E. coli., yeast, insect cells, and mammalian cells.




Fig. 1-1 Schematic diagram of ligand of IMAC medium.

### ■ IDA

Ni IDA Beads and Ni IDA Beads 6FF are coupled with iminodiacetic acid (IDA), forming a stable planar quadrilateral structure after chelating with precharged nickel ions; this structure can be easily disrupted by competitive small molecules. The Ni-IDA structure provides multiple binding sites for histidine-tagged proteins. Ni IDA Beads 6FF offers advantages in terms of higher binding capacity and cost-effectiveness, along with excellent flow properties that facilitate scaling up.

### 📦 Product Information

| Products         | Cat.No.  | Size        | Applications   |
|------------------|----------|-------------|--|
| Ni IDA Beads     | SA003005 | 5ml         | Ni IDA Beads and 6FF are cost-effective products for purifying polyhistidine-tagged proteins and are recommended for the purification of soluble protein extracts. |
|                  | SA003025 | 25ml        |  |
|                  | SA003100 | 100ml       |  |
|                  | SA003500 | 500ml       |  |
|                  | SA00301L | 1L          |  |
|                  | SA00310L | 10L         |  |
| Ni IDA Beads 6FF | SA052005 | 5ml         |  |
|                  | SA052025 | 25ml        |  |
|                  | SA052100 | 100ml       |  |
|                  | SA052500 | 500ml       |  |
|                  | SA05201L | 1L          |  |
|                  | SA05210L | 10L         |  |
| HisCap IDA 6FF   | SA052C11 | 1×1ml       |  |
|                  | SA052C15 | 1×5ml       |  |
|                  | SA052C51 | 5×1ml       |  |
|                  | SA052C55 | 5×5ml       |  |
|                  | SA052CS  | 3×1ml+1×5ml |  |

NTA

Ni NTA Beads, Ni NTA Beads 6FF, Ni NTA Beads HF utilize nitrilotriacetic acid (NTA) as the ligand. Nickel ions (Ni<sup>2+</sup>) are chelated to this group, resulting in a stable octahedral structure that protects the nickel ions from competitive small molecule interference. The Ni-NTA structure is compatible with certain concentrations of reducing agents, denaturing agents, detergents, and other additives.

Product Information

| Products                        | Cat.No          | Size        | Applications   |  |
|---------------------------------|-----------------|-------------|--|--|
| Ni NTA Beads                    | SA004005        | 5ml         | Ni NTA Beads can be used in the purification of polyhistidine-tagged inclusion bodies and soluble proteins. It can withstand more severe chemical reagents and has a wider range of applications.                                  |  |
|                                 | SA004025        | 25ml        |  |  |
|                                 | SA004100        | 100ml       |  |  |
|                                 | SA004500        | 500ml       |  |  |
|                                 | SA00401L        | 1L          |  |  |
|                                 | SA00410L        | 10L         |  |  |
| Ni NTA Beads Gravity Column     | SA004GC01       | 1ml         | It is cost-effective products for purifying His-tagged proteins, offering simple operation and high purification efficiency.   |  |
| HisPur Ni NTA Kit               | SA004K03        | 3T          |  |  |
| HisPur Ni NTA Buffer Kit        | SA004KB03       | 3T          | Ni NTA Beads 6FF is compatible with both native and denaturing conditions, making it suitable for purifying polyhistidine-tagged inclusion bodies and soluble proteins. Its high flow properties make it excellent for scaling-up. |  |
| Ni NTA Beads 6FF                | SA005005        | 5ml         |  |  |
|                                 | SA005025        | 25ml        |  |  |
|                                 | SA005100        | 100ml       |  |  |
|                                 | SA005500        | 500ml       |  |  |
|                                 | SA00501L        | 1L          |  |  |
|                                 | SA00510L        | 10L         |  |  |
| Ni NTA Beads 6FF Gravity Column | SA005GC01       | 1ml         |  | It can be used in the purification of polyhistidine-tagged inclusion bodies and soluble proteins. Highly rigid agarose microspheres withstand elevated flow rates and deliver enhanced chemical stability. |
| HisCap 6FF                      | SA005C11        | 1x1ml       |  |  |
| HisCap 6FF                      | SA005C51        | 5x1ml       |  |  |
|                                 | SA005C15        | 1x5ml       |  |  |
|                                 | SA005C55        | 5x5ml       |  |  |
|                                 | SA005CS         | 3x1ml+1x5ml |  |  |
|                                 | HiSelect Ni 6FF | SA005C47    | 4.7ml, 1pcs  |  |
| HiPur Ni 6FF                    | SA005C20        | 20ml, 1pcs  | It can be used in the purification of polyhistidine-tagged inclusion bodies and soluble proteins. Highly rigid agarose microspheres withstand elevated flow rates and deliver enhanced chemical stability.                         |  |
| Ni NTA Beads HF                 | SA101005        | 5ml         |  |  |
|                                 | SA101025        | 25ml        |  |  |
|                                 | SA101100        | 100ml       |  |  |
|                                 | SA101500        | 500ml       |  |  |
|                                 | SA10101L        | 1L          |  |  |
|                                 | SA10110L        | 10L         |  |  |
| Ni NTA Beads HF Gravity Column  | SA100GC01       | 1ml         |  | It can be used in the purification of polyhistidine-tagged inclusion bodies and soluble proteins. Highly rigid agarose microspheres withstand elevated flow rates and deliver enhanced chemical stability. |
| HisCap NTA Beads HF             | SA101C11        | 1x1ml       |  |  |
| HisCap NTA Beads HF             | SA101C51        | 5x1ml       |  |  |
|                                 | SA101C15        | 1x5ml       |  |  |
|                                 | SA101C55        | 5x5ml       |  |  |
|                                 | HiSelect Ni HF  | SA101C47    | 4.7ml, 1pcs  |  |
| HiPur Ni HF                     | SA101C20        | 20ml, 1pcs  |  |  |

Smart

Ni Smart Beads, Ni Smart Beads 6FF and Ni Smart Beads HF are novel immobilized metal ion affinity chromatography (IMAC) media precharged with nickel ions, specifically designed for the capture and purification of histidine-tagged proteins secreted into eukaryotic cell culture supernatants. The robust nickel ion binding provides exceptional resistance to EDTA and reducing agents such as DTT, allowing for the direct loading of large sample volumes without the need to remove agents that typically strip nickel ions from conventional IMAC media. Ni Smart Beads 6FF are stable in all buffers commonly used in IMAC; refer to Table 1-2 for information on reagent tolerance.

Example: Ni Smart Beads 6FF

Ni Smart Beads 6FF can be directly utilized for the capture and purification of 6xHis-tagged proteins expressed by mammalian cells, without the need to change the solution or concentrate the cellular supernatant.




Fig.1-2 Purification of 6xHis tagged protein expressed by mammalian cells.

His-tag usually consists of 6-10 consecutive histidine residues. The structure of the tag (including its position, sequence and length) can affect various aspects of protein production, such as expression levels, three-dimensional structure, crystallization, and binding capacity.




Fig.1-3 Purification Results of Proteins with Different His-Tags

Fig.1-4 SDS-PAGE Analysis of 6xHis-Tagged Protein Purification

Product Information

| Products                          | Cat.No    | Size        | Applications   |
|-----------------------------------|-----------|-------------|--|
| Ni Smart Beads                    | SA035005  | 5ml         | Ni Smart Beads and 6FF offer exceptional resistance to a wide range of chemicals, making them suitable for protein environments with high concentrations of reducing and chelating agents. They are especially effective in purification systems for proteins expressed in eukaryotic systems.   |
|                                   | SA035025  | 25ml        |  |
|                                   | SA035100  | 100ml       |  |
|                                   | SA035500  | 500ml       |  |
|                                   | SA03501L  | 1L          |  |
|                                   | SA03510L  | 10L         |  |
| Ni Smart Beads Gravity Column     | SA035GC01 | 1ml         |  |
|                                   | SA035GC05 | 5ml         |  |
| Ni Smart Beads 6FF                | SA036005  | 5ml         |  |
|                                   | SA036025  | 25ml        |  |
|                                   | SA036100  | 100ml       |  |
|                                   | SA036500  | 500ml       |  |
|                                   | SA03601L  | 1L          |  |
|                                   | SA03610L  | 10L         |  |
| Ni Smart Beads 6FF Gravity Column | SA036GC01 | 1ml         |  |
|                                   | SA036GC05 | 5ml         |  |
| HisCap Smart 6FF                  | SA036C11  | 1×1ml       |  |
|                                   | SA036C51  | 5×1ml       |  |
|                                   | SA036C15  | 1×5ml       |  |
|                                   | SA036C55  | 5×5ml       |  |
|                                   | SA036CS   | 3×1ml+1×5ml |  |
| HiSelect Ni Smart                 | SA036C47  | 4.7ml,1pcs  |  |
| HiPur Ni Smart                    | SA036C20  | 20ml,1pcs   |  |
| Ni Smart Beads HF                 | SA110005  | 5ml         | It offer exceptional resistance to a wide range of chemicals, making them suitable for protein environments with high concentrations of reducing and chelating agents. They are especially effective in purification systems for proteins expressed in eukaryotic systems.Highly rigid agarose microspheres withstand elevated flow rates and deliver enhanced chemical stability. |
|                                   | SA110025  | 25ml        |  |
|                                   | SA110100  | 100ml       |  |
|                                   | SA110500  | 500ml       |  |
|                                   | SA11001L  | 1L          |  |
| SA11010L                          | 10L       |             |  |
| Ni Smart Beads HF Gravity Column  | SA110GC01 | 1ml         |  |
|                                   | SA110GC05 | 5ml         |  |
| HisCap Smart Beads HF             | SA110C11  | 1×1ml       |  |
|                                   | SA110C51  | 5×1ml       |  |
|                                   | SA110C15  | 1×5ml       |  |
|                                   | SA110C55  | 5×5ml       |  |
| HiSelect Ni HF                    | SA110C47  | 4.7ml,1pcs  |  |
| HiPur Ni Smart HF                 | SA110C20  | 20ml,1pcs   |  |
| HiPur Ni Smart Kit                | SA110K    | 5T          |  |

Other metal ion chelating affinity purification resin

Chelating Beads 6FF and IMAC Beads 6FF are highly cross-linked 6% agarose beads functionalized with IDA and NTA ligands, enabling them to chelate various metal ions, such as Cu<sup>2+</sup>, Zn<sup>2+</sup>, Co<sup>2+</sup>, Ni<sup>2+</sup>, and Fe<sup>3+</sup>. When choosing the desired metal ion, it is important to consider the structural requirements for optimal metal chelate-protein binding. The binding affinity to polyhistidine tags generally follows the order: Cu<sup>2+</sup> > Ni<sup>2+</sup> > Zn<sup>2+</sup> > Co<sup>2+</sup>. Typically, Cu<sup>2+</sup> binds strongly to histidine-tagged proteins and is effective for purifying proteins with weaker interactions, while Co<sup>2+</sup> provides greater specificity for histidine tags, yielding higher-purity samples.

Product Information

| Products                        | Cat.No      | Size       | Applications  |
|---------------------------------|-------------|------------|---|
| Cu IDA Beads                    | SA041005    | 5ml        | Cu IDA Beads are beads with an average particle size of 90 µm coupled with IDA ligands. The medium precharged with Cu <sup>2+</sup> ions.   |
|                                 | SA041025    | 25ml       |   |
|                                 | SA041100    | 100ml      |   |
|                                 | SA041500    | 500ml      |   |
|                                 | SA04101L    | 1L         |   |
|                                 | SA04110L    | 10L        |   |
| Cu IDA Beads Gravity Column     | SA041GC01   | 1ml        |   |
|                                 | SA041GC05   | 5ml        |   |
| Zn IDA Beads                    | SA080005    | 5ml        | Zn IDA Beads are beads with an average particle size of 90 µm coupled with IDA ligands. The medium precharged with Zn <sup>2+</sup> .   |
|                                 | SA080025    | 25ml       |   |
|                                 | SA080100    | 100ml      |   |
|                                 | SA080500    | 500ml      |   |
|                                 | SA08001L    | 1L         |   |
|                                 | SA08010L    | 10L        |   |
| Co NTA Beads                    | SA037005    | 5ml        | Co NTA Beads are beads with an average particle size of 90 µm coupled with NTA ligands. The medium is precharged with Co <sup>2+</sup> ions.  |
|                                 | SA037025    | 25ml       |   |
|                                 | SA037100    | 100ml      |   |
|                                 | SA037500    | 500ml      |   |
|                                 | SA03701L    | 1L         |   |
| Co NTA Beads Gravity Column     | SA03710L    | 10L        |   |
|                                 | SA037GC01   | 1ml        |   |
| Co NTA Beads Gravity Column     | SA037GC05   | 5ml        |   |
|                                 | SA038005    | 5ml        | Co NTA Beads 6FF is compatible with both native and denaturing conditions, making them suitable for the purification of polyhistidine-tagged inclusion bodies as well as soluble proteins. Its high flow properties make it excellent for scaling-up. |
| SA038025                        | 25ml        |            |   |
| SA038100                        | 100ml       |            |   |
| SA038500                        | 500ml       |            |   |
| SA03801L                        | 1L          |            |   |
| Co NTA Beads 6FF                | SA03810L    | 10L        |   |
|                                 | SA038GC01   | 1ml        |   |
| Co NTA Beads 6FF Gravity Column | SA038GC05   | 5ml        |   |
|                                 | SA038C11    | 1×1ml      | The columns have the standard interface that can be adapted to the liquid chromatography systems, such as Smart GO™ and Smart PUR™.   |
| SA038C51                        | 5×1ml       |            |   |
| SA038C15                        | 1×5ml       |            |   |
| SA038C55                        | 5×5ml       |            |   |
| SA038CS                         | 3×1ml+1×5ml |            |   |
| HiSelect Co 6FF                 | SA038C47    | 4.7ml,1pcs |   |
| HiPur Co 6FF                    | SA038C20    | 20ml,1pcs  |   |
| IMAC Beads 6FF                  | SA050005    | 5ml        | IMAC Beads 6FF is the medium without chelating metal ion and its ligand is NTA. It is compatible with a certain concentration of reducing agents, denaturing agents, detergents and other additives.  |
|                                 | SA050025    | 25ml       |   |
|                                 | SA050100    | 100ml      |   |
|                                 | SA050500    | 500ml      |   |
|                                 | SA05001L    | 1L         |   |
|                                 | SA05010L    | 10L        |   |
| Chelating Beads 6FF             | SA051005    | 5ml        | Chelating Beads 6FF is the medium without chelating metal ion and its ligand is IDA.  |
|                                 | SA051025    | 25ml       |   |
|                                 | SA051100    | 100ml      |   |
|                                 | SA051500    | 500ml      |   |
|                                 | SA05101L    | 1L         |   |
|                                 | SA05110L    | 10L        |   |

## GST fusion protein affinity purification resin

Glutathione Beads and Glutathione Beads 4FF are affinity chromatography media designed for the purification of glutathione S-transferase (GST), glutathione-dependent proteins, and GST fusion proteins expressed in various systems. The reduced glutathione ligand is coupled to agarose beads via a 12-atom linker, which is optimized to achieve a high binding capacity for GST fusion proteins. Glutathione Beads 4FF can be utilized in fast protein liquid chromatography (FPLC) and is easy to pack and use, with high flow properties that facilitate scale-up. The GST Pur Glutathione Kit is designed for the one-step purification of GST fusion proteins and consists of a prepacked column and ready-to-use buffer. This kit eliminates the need for time-consuming buffer preparation, thus promoting rapid and convenient purification of the target protein from complex samples.

### Example: Glutathione Beads




Fig. 1-5 SDS-PAGE analysis of purification of GST fusion protein.

### Product Information

| Products                             | Cat.No.   | Size        | Applications  |
|--------------------------------------|-----------|-------------|---|
| Glutathione Beads                    | SA008005  | 5ml         | Glutathione Beads offer a high binding capacity for GST fusion proteins and are recommended for laboratory-scale                                  |
|                                      | SA008025  | 25ml        |   |
|                                      | SA008100  | 100ml       |   |
|                                      | SA008500  | 500ml       |   |
|                                      | SA00801L  | 1L          |   |
| Glutathione Beads Gravity Column     | SA00810L  | 10L         | The GST fusion protein purification kit offers simple operation and high purification efficiency.   |
|                                      | SA008GC01 | 1ml         |   |
|                                      | SA008GC05 | 5ml         |   |
| GSTPur Glutathione Kit               | SA008K03  | 3T          | The columns have the standard interface that can be adapted to the liquid chromatography systems, such as Smart GO™ and Smart PUR™.               |
|                                      | SA008K05  | 5T          |   |
| GSTPur Glutathione Buffer Kit        | SA008KB03 | 3T          | The columns have the standard interface that can be adapted to the liquid chromatography systems, such as Smart GO™ and Smart PUR™.               |
|                                      | SA008KB05 | 5T          |   |
| Glutathione Beads 4FF                | SA010005  | 5ml         | Glutathione Beads 4FF can be used in purification systems and is easy to pack and use. Its high flow properties make it excellent for scaling-up. |
|                                      | SA010025  | 25ml        |   |
|                                      | SA010100  | 100ml       |   |
|                                      | SA010500  | 500ml       |   |
|                                      | SA01001L  | 1L          |   |
| Glutathione Beads 4FF Gravity Column | SA01010L  | 10L         | The columns have the standard interface that can be adapted to the liquid chromatography systems, such as Smart GO™ and Smart PUR™.               |
|                                      | SA010GC01 | 1ml         |   |
|                                      | SA010GC05 | 5ml         |   |
| GSTCap 4FF                           | SA010C11  | 1×1ml       | The columns have the standard interface that can be adapted to the liquid chromatography systems, such as Smart GO™ and Smart PUR™.               |
|                                      | SA010C51  | 5×1ml       |   |
|                                      | SA010C15  | 1×5ml       |   |
|                                      | SA010C55  | 5×5ml       |   |
|                                      | SA010CS   | 3×1ml+1×5ml |   |
| HiSelect Glutathione 4FF             | SA010C47  | 4.7ml, 1pcs |   |
| HiPur Glutathione 4FF                | SA010C20  | 20ml, 1pcs  |   |

## MBP fusion protein affinity purification resin

Maltose-binding protein (MBP) is a 42 kDa tag that enhances the soluble expression of target proteins, promoting correct folding and increasing expression levels. Dextrin Beads 6FF serves as an affinity chromatography medium for purifying MBP fusion proteins, enabling one-step purification. Bound fusion proteins can be gently eluted with 10 mM maltose to maintain the target protein's activity. MBP can be removed from the fusion protein using site-specific protease.

### Example: Dextrin Beads 6FF




Fig. 1-6 SDS-PAGE analysis of purification of MBP fusion protein.

### Product Information

| Products                         | Cat.No.   | Size        | Applications   |
|----------------------------------|-----------|-------------|--|
| Dextrin Beads                    | SA077005  | 5ml         | Dextrin Beads are suitable for purifying MBP and MBP fusion proteins, offering high purity and strong binding capacity for the target protein. |
|                                  | SA077025  | 25ml        |  |
|                                  | SA077100  | 100ml       |  |
|                                  | SA077500  | 500ml       |  |
|                                  | SA07701L  | 1L          |  |
| Dextrin Beads Gravity Column     | SA07710L  | 10L         | Its high flow properties make it excellent for scaling-up.   |
|                                  | SA077GC01 | 1ml         |  |
|                                  | SA077GC05 | 5ml         |  |
| Dextrin Beads 6FF                | SA026005  | 5ml         | The columns have the standard interface that can be adapted to the liquid chromatography systems.  |
|                                  | SA026025  | 25ml        |  |
|                                  | SA026100  | 100ml       |  |
|                                  | SA026500  | 500ml       |  |
|                                  | SA02601L  | 1L          |  |
| Dextrin Beads 6FF Gravity Column | SA02610L  | 10L         | The columns have the standard interface that can be adapted to the liquid chromatography systems.  |
|                                  | SA026GC01 | 1ml         |  |
|                                  | SA026GC05 | 5ml         |  |
| PreCap Dextrin                   | SA026C11  | 1×1ml       | The columns have the standard interface that can be adapted to the liquid chromatography systems.  |
|                                  | SA026C51  | 5×1ml       |  |
|                                  | SA026C15  | 1×5ml       |  |
|                                  | SA026C55  | 5×5ml       |  |
|                                  | SA026CS   | 3×1ml+1×5ml |  |
| HiSelect Dextrin 6FF HiPur       | SA026C47  | 4.7ml, 1pcs |  |
| Dextrin 6FF                      | SA026C20  | 20ml, 1pcs  |  |

### Biotin-tagged protein affinity purification resin

Streptavidin Beads 6FF is an affinity chromatography medium used for biotin-streptavidin interactions. The interaction between streptavidin and biotin is extremely strong, requiring denaturing conditions for elution, which can damage both the ligand and the sample. However, the weaker interaction between 2-iminobiotin and streptavidin allows for elution at pH 4.0, preserving biomolecular activity. The crosslinking of the base matrix has been optimized to ensure high physical and chemical stability, making it suitable for cost-effective, large-scale applications.

#### Example: Streptavidin Beads 6FF




Fig. 1-7 SDS-PAGE analysis of purification of iminobiotinylated protein.

### Strep-tag II fusion protein affinity purification resin

Streptavidin ST Beads 4FF is a chromatography medium designed for one-step purification of Strep-tag II fusion proteins from various expression systems. The Strep-tag II is an eight-amino acid peptide (Trp-Ser-His-Pro-Gln-Phe-Glu-Lys), which don't need to be removed due to their minimal effect on fusion proteins. The Streptactin ligand, immobilized on highly cross-linked 4% agarose, is a specially recombinant protein with a binding affinity for Strep-tag II that exceeds that of streptavidin. Purification occurs under physiological conditions, and mild elution preserves the target protein's activity.

#### Example: STarm Beads 4FF




Fig. 1-8 Effect of Biotin Concentration on Elution Efficiency of STarm Beads 4FF

### Product Information

| Products                              | Cat.No.   | Size        | Applications   |
|---------------------------------------|-----------|-------------|--|
| Streptavidin Beads 6FF                | SA021005  | 5ml         | Streptavidin Beads 6FF can be used for the immobilization or detection of biotinylated nucleotides, peptides, proteins, and cells.                 |
|                                       | SA021025  | 25ml        |  |
|                                       | SA021100  | 100ml       |  |
|                                       | SA021500  | 500ml       |  |
|                                       | SA02101L  | 1L          |  |
|                                       | SA02110L  | 10L         |  |
| Streptavidin Beads 6FF Gravity Column | SA021GC01 | 1ml         |  |
|                                       | SA021GC05 | 5ml         |  |
| PreCap Streptavidin                   | SA021C11  | 1×1ml       | PreCap Streptavidin columns have the standard interface and can be adapted to the liquid chromatography systems, such as Smart GO™ and Smart PUR™. |
|                                       | SA021C51  | 5×1ml       |  |
|                                       | SA021C15  | 1×5ml       |  |
|                                       | SA021C55  | 5×5ml       |  |
|                                       | SA021CS   | 3×1ml+1×5ml |  |
| HiSelect Streptavidin 6FF             | SA021C47  | 4.7ml,1pcs  |  |
| HiPur Streptavidin 6FF                | SA021C20  | 20ml,1pcs   |  |

### Product Information

| Products                     | Cat.No.  | Size        | Applications  |
|------------------------------|----------|-------------|---|
| Streptavidin ST Beads 4FF    | SA053005 | 5ml         | Streptactin Beads 4FF is a chromatography medium for one-step purifying Strep-tag II fusion proteins.   |
|                              | SA053025 | 25ml        |   |
|                              | SA053100 | 100ml       |   |
|                              | SA053500 | 500ml       |   |
|                              | SA05301L | 1L          |   |
| PreCap Streptavidin ST       | SA05310L | 10L         |   |
|                              | SA053C11 | 1×1ml       |   |
|                              | SA053C51 | 5×1ml       |   |
|                              | SA053C15 | 1×5ml       |   |
|                              | SA053C55 | 5×5ml       |   |
|                              | SA053CS  | 3×1ml+1×5ml |   |
| HiSelect Streptavidin ST 4FF | SA053C47 | 4.7ml,1pcs  | STarm Streptactin Beads 4FF is uniquely designed to purify Strep-tag II or Twin Strep-Tag II tagged proteins. These beads are widely used for high load and gentle regeneration conditions. It can be reused for many times and has stronger tolerance. |
| HiPur Streptavidin ST 4FF    | SA053C20 | 20ml,1pcs   |   |
| STarm Beads 4FF              | SA092005 | 5ml         |   |
|                              | SA092025 | 25ml        |   |
|                              | SA092100 | 100ml       |   |
|                              | SA092500 | 500ml       |   |
|                              | SA09201L | 1L          |   |
| PreCap STarm                 | SA09210L | 10L         |   |
|                              | SA092C11 | 1×1ml       |   |
|                              | SA092C51 | 5×1ml       |   |
|                              | SA092C15 | 1×5ml       |   |
|                              | SA092C55 | 5×5ml       |   |
|                              | SA092CS  | 3×1ml+1×5ml |   |
| HiSelect STarm 4FF HiPur     | SA092C47 | 4.7ml,1pcs  |   |
| STarm 4FF                    | SA092C20 | 20ml,1pcs   |   |

## DYKDDDDK-tagged protein affinity purification resin

Anti-DYKDDDDK Affinity Beads is used for the one-step purification of DYKDDDDK (Flag)-tagged proteins expressed in prokaryotic, yeast, or mammalian cells. The Flag sequence consists of eight hydrophilic amino acids and can be removed by enterokinase. The medium is capable of binding the Flag tag at the N-terminal and C-terminal of the fusion protein. The ligand of the medium is an anti-DYKDDDDK antibody that is covalently coupled to 4% agarose beads. Anti-DYKDDDDK Affinity Beads exhibit high specificity for the fusion protein and can be used for both purification and immunoprecipitation (IP) of Flag-tagged proteins.

### Example: Anti-DYKDDDDK HC Affinity Beads

#### 1) Prokaryotic Protein Purification :




Fig. 1-9. SDS-PAGE: Prokaryotic Samples Purified with Anti-DYKDDDDK HC Affinity Beads

#### 2) Eukaryotic Protein Purification:




Fig. 1-10. SDS-PAGE: Eukaryotic Samples Purified with Anti-DYKDDDDK HC Affinity Beads

#### Product Information

| Products                     | Cat.No.  | Size  | Applications   |
|------------------------------|----------|-------|--|
| Anti-DYKDDDDK Affinity Beads | SA042001 | 1ml   | One step purification and detection of the DYKDDDDK-tagged protein expressed in prokaryotic, yeast or mammalian cells. |
|                              | SA042005 | 5ml   |  |
|                              | SA042025 | 25ml  |  |
|                              | SA042100 | 100ml |  |
|                              | SA042500 | 500ml |  |
|                              | SA04201L | 1L    |  |

| Products                        | Cat.No.  | Size  | Applications  |
|---------------------------------|----------|-------|---|
| Anti-DYKDDDDK SI Affinity Beads | SA109001 | 1ml   | One step purification and detection of the DYKDDDDK-tagged protein expressed in prokaryotic, yeast or mammalian cells.                                  |
|                                 | SA109005 | 5ml   |   |
|                                 | SA109025 | 25ml  |   |
|                                 | SA109100 | 100ml |   |
|                                 | SA109500 | 500ml |   |
| Anti-DYKDDDDK HC Affinity Beads | SA10901L | 1L    | Delivers excellent purification performance for both prokaryotic and eukaryotic samples, while tolerating relatively high concentrations of detergents. |
|                                 | SA145001 | 1ml   |   |
|                                 | SA145005 | 5ml   |   |
|                                 | SA145025 | 25ml  |   |
|                                 | SA145100 | 100ml |   |
|                                 | SA145500 | 500ml |   |
|                                 | SA14501L | 1L    |   |

## HA-tagged protein affinity purification resin

The HA tag is a peptide consisting of the 98th to 106th amino acids (YPYDVPDYA) of human influenza hemagglutinin and has minimal impact on the spatial structure of the target protein. This peptide can be fused to either the N-terminal or C-terminal of the protein, making it commonly used in the expression of recombinant proteins. Anti-HA Affinity Beads are based on 4% agarose beads and demonstrate high specificity for binding HA-tagged proteins. These beads can be used for the purification and immunoprecipitation (IP) of HA-tagged proteins.

#### Product Information

| Products               | Cat.No.  | Size  | Applications  |
|------------------------|----------|-------|---|
| Anti-HA Affinity Beads | SA068001 | 1ml   | It is used for purification and detection of HA-tagged protein. |
|                        | SA068005 | 5ml   |   |
|                        | SA068025 | 25ml  |   |
|                        | SA068100 | 100ml |   |
|                        | SA068500 | 500ml |   |
|                        | SA06801L | 1L    |   |

## c-Myc-tagged protein affinity purification resin

Anti-c-Myc Affinity Beads are used for the detection and purification of c-Myc-tagged proteins expressed in both prokaryotic and eukaryotic cells. The c-Myc peptide is encoded by the c-Myc gene located on human chromosome 8q24, corresponding to the 410th to 419th amino acids (EQKLISEEDL) of human p62. These affinity beads can recognize c-Myc-tagged proteins at the C-terminal, N-terminal, or internal sites. They are suitable for applications such as Western blot hybridization, immunoprecipitation, and flow cytometry.

#### Product Information

| Products                  | Cat.No.  | Size  | Applications   |
|---------------------------|----------|-------|--|
| Anti-c-Myc Affinity Beads | SA065001 | 1ml   | Anti-c-Myc Affinity Beads can be used to detect and purify the c-myc tagged proteins expressed by prokaryotic and eukaryotic cell. |
|                           | SA065005 | 5ml   |  |
|                           | SA065025 | 25ml  |  |
|                           | SA065100 | 100ml |  |
|                           | SA065500 | 500ml |  |
|                           | SA06501L | 1L    |  |

## GFP-tagged protein affinity purification resin

Green fluorescent protein (GFP) and its mutant, enhanced green fluorescent protein (EGFP), are widely used to detect the expression efficiency and distribution of target genes and proteins. Both GFP and GFP-tagged proteins can be detected spontaneously through fluorescence, requiring no special procedures, and are unaffected by other substances. Anti-GFP Affinity Beads 4FF can be used to detect and purify GFP, EGFP, and their fusion proteins without binding to blue fluorescent protein (BFP)-tagged proteins.




Fig.1-11 WB analysis of GFP-tagged protein.




Fig.1-12 SDS-PAGE analysis of purification of GFP-tagged protein.

### Product Information

| Products                    | Cat.No.  | Size  | Applications   |
|-----------------------------|----------|-------|--|
| Anti-GFP Affinity Beads 4FF | SA070001 | 1ml   | Anti-GFP Affinity Beads 4FF can be used to detect and purify GFP, EGFP and their fusion protein without binding to BFP fusion protein. |
|                             | SA070005 | 5ml   |  |
|                             | SA070025 | 25ml  |  |
|                             | SA070100 | 100ml |  |
|                             | SA070500 | 500ml |  |
|                             | SA07001L | 1L    |  |

## RFP-tagged protein affinity purification resin

Red fluorescent protein (RFP), a molecular marker isolated from anemones, is commonly used in studies of protein localization, gene expression, and protein interactions. Anti-RFP Affinity Beads 4FF, based on cross-linked agarose beads, can be used for immunoprecipitation, co-immunoprecipitation, and the purification of natural RFP, RFP mutants, and their fusion proteins.

### Product Information

| Products                    | Cat.No.  | Size  | Applications  |
|-----------------------------|----------|-------|---|
| Anti-RFP Affinity Beads 4FF | SA072001 | 1ml   | Anti-RFP Affinity Beads 4FF can be used for detection and purification of RFP and its fusion protein. |
|                             | SA072005 | 5ml   |   |
|                             | SA072025 | 25ml  |   |
|                             | SA072100 | 100ml |   |
|                             | SA072500 | 500ml |   |
|                             | SA07201L | 1L    |   |

## YFP-tagged protein affinity purification resin

Yellow fluorescent protein (YFP) is a genetic variant of green fluorescent protein (GFP). The YFP-tagged protein fluoresces spontaneously, allowing for the localization of the target protein within the cell, and with minimal interference from other substances. This product is applicable for IP, co-IP, and affinity chromatography of natural YFP and its fusion proteins. Additionally, Anti-YFP Affinity Beads 4FF can detect and purify GFP, EGFP, and their fusion proteins, while not binding to BFP or RFP-tagged proteins.

### Product Information

| Products                    | Cat.No.  | Size  | Applications   |
|-----------------------------|----------|-------|--|
| Anti-YFP Affinity Beads 4FF | SA075001 | 1ml   | Anti-YFP Affinity Beads 4FF can be used for immunoprecipitation and affinity chromatography of YFP and its fusion protein. |
|                             | SA075005 | 5ml   |  |
|                             | SA075025 | 25ml  |  |
|                             | SA075100 | 100ml |  |
|                             | SA075500 | 500ml |  |
|                             | SA07501L | 1L    |  |

## Protein A fusion protein affinity purification resin

Rabbit IgG Beads 4FF is a one-step affinity chromatography for the purification of protein A fusion proteins. Protein A, a cell wall protein derived from *Staphylococcus aureus*, specifically binds to the Fc fragment of mammalian IgG. Consequently, Rabbit IgG is coupled to cross-linked 4% agarose beads, enabling effective purification and detection of protein A fusion proteins.

### Product Information

| Products         | Cat.No.  | Size       | Applications   |
|------------------|----------|------------|--|
| Rabbit IgG Beads | SA082005 | 5ml        | Rabbit IgG Beads 4FF can be used to purify and detect the protein A fusion proteins expressed by prokaryotic and eukaryotic cells. |
|                  | SA082025 | 25ml       |  |
|                  | SA082100 | 100ml      |  |
|                  | SA082500 | 500ml      |  |
|                  | SA08201L | 1L         |  |
| HiSelect IgG 4FF | SA028C47 | 4.7ml,1pcs |  |
| HiPur IgG 4FF    | SA028C20 | 20ml,1pcs  |  |

# Antibody Affinity Chromatography

The ligands used in the purification of general antibodies mainly include Protein A and Protein G. These ligands are commonly used for the purification of both polyclonal and monoclonal antibodies.




Fig. 1-13 the site of the antibody bind to the protein.

Table 3-1 The binding ability of protein A, protein G and protein A/G to different antibodies.

| Species         | Subtype   | Protein A | Protein G | Protein A/G |
|-----------------|-----------|-----------|-----------|-------------|
| Human           | IgA       | variable  | -         | ++          |
|                 | IgD       | -         | -         | -           |
|                 | IgE       | -         | -         | -           |
|                 | IgG1      | ++++      | ++++      | ++++        |
|                 | IgG2      | ++++      | ++++      | ++++        |
|                 | IgG3      | -         | ++++      | ++++        |
|                 | IgG4      | ++++      | ++++      | ++++        |
| Avian egg yolk  | IgM       | variable  | -         | ++          |
|                 | IgY       | -         | -         | -           |
| Cow             | -         | ++        | ++++      | ++++        |
| Dog             | -         | ++++      | ++        | ++++        |
| Goat            | -         | -         | ++++      | ++++        |
| Guinea          | IgG1      | ++++      | ++        | ++++        |
|                 | IgG2      | ++++      | ++        | ++++        |
| pig             | -         | +         | ++        | -           |
| Hamster Horse   | Total IgG | ++        | ++++      | ++++        |
| Koala           | -         | -         | +         | -           |
| Llama           | -         | -         | +         | -           |
| Monkey (rhesus) | IgG1      | ++++      | ++++      | ++++        |
|                 | IgG2a     | +         | ++++      | ++          |
| Mouse           | IgG2b     | ++++      | ++++      | ++++        |
|                 | IgG3      | +++       | +++       | +++         |
|                 | IgM       | ++        | ++        | ++          |
|                 | -         | variable  | -         | -           |
| Pig             | Total     | +++       | +++       | ++++        |
| Rabbit          | IgG       | ++++      | +++       | ++++        |
|                 | IgG1      | -         | +         | ++          |
| Rat             | IgG2a     | -         | ++++      | ++++        |
|                 | IgG2b     | -         | ++        | ++          |
|                 | IgG3      | +         | ++        | ++          |
| Sheep           | Total     | +/-       | ++        | ++          |

++++=Strong binding; ++ = medium binding; - = weak or no binding

## Protein A related products

### Universal antibody purification resin

Protein A, a cell wall protein isolated from *Staphylococcus aureus*, mainly binds to mammalian IgG via the Fc fragment but does not bind to human IgM, IgD, or IgA. Natural Protein A contains five IgG binding regions along with several unknown functional regions, while recombinant Protein A retains only the five IgG binding regions, which has low nonspecific adsorption.

The AbPur rProtein A Kit is tailored for the one-step purification of antibodies, which contains one pre-packed column and ready-to-use buffers.

### Product Information

| Products                        | Cat.No.   | Size       | Applications   |
|---------------------------------|-----------|------------|--|
| rProtein A Beads                | SA012005  | 5ml        | The kit can be used to purify both monoclonal and polyclonal antibodies, and the purity can reach 95% in one step. It is recommended for use at the laboratory scale.                |
|                                 | SA012025  | 25ml       |  |
|                                 | SA012100  | 100ml      |  |
|                                 | SA012500  | 500ml      |  |
|                                 | SA01201L  | 1L         |  |
| rProtein A Beads Gravity Column | SA012GC01 | 1ml        |  |
|                                 | SA012GC05 | 5ml        |  |
| AbPur rProtein A Kit            | SA012K03  | 3T         |  |
| AbPur rProtein A Buffer Kit     | SA012KB03 | 3T         |  |
|                                 | SA015005  | 5ml        |  |
| rProtein A Beads 4FF            | SA015025  | 25ml       | rProtein A Beads 4FF can be used in FPLC and is easy to pack and use. Its high flow properties make it excellent for scaling-up.   |
|                                 | SA015100  | 100ml      |  |
|                                 | SA015500  | 500ml      |  |
|                                 | SA01501L  | 1L         |  |
|                                 | SA01510L  | 10L        |  |
| rProtein A Beads                | SA015GC01 | 1ml        |  |
|                                 | SA015GC05 | 5ml        |  |
| AbCap A 4FF                     | SA015C11  | 1×1ml      | The columns have the standard interface that can be adapted to the liquid chromatography systems, such as Smart GO™ and Smart PUR™.  |
|                                 | SA015C51  | 5×1ml      |  |
|                                 | SA015C15  | 1×5ml      |  |
|                                 | SA015C55  | 5×5ml      |  |
| HiSelect A 4FF                  | SA015C47  | 4.7ml,1pcs |  |
|                                 | SA015C20  | 20ml,1pcs  |  |
| Smac A-ME                       | SA123005  | 5ml        | Enables elution of Fc-containing molecules under moderate pH conditions.   |
|                                 | SA123025  | 25ml       |  |
|                                 | SA123100  | 100ml      |  |
|                                 | SA123500  | 500ml      |  |
|                                 | SA12301L  | 1L         |  |
| PreCap Smac A-ME                | SA12310L  | 10L        | The columns have the standard interface that can be adapted to the liquid chromatography systems, such as Smart GO™ and Smart PUR™.  |
|                                 | SA123C11  | 1×1ml      |  |
|                                 | SA123C51  | 5×1ml      |  |
|                                 | SA123C15  | 1×5ml      |  |
| Smac VH3 Affinity Beads         | SA123C55  | 5×5ml      |  |
|                                 | SA122005  | 5ml        |  |
|                                 | SA122025  | 25ml       |  |
|                                 | SA122100  | 100ml      |  |
|                                 | SA122500  | 500ml      |  |
| PreCap Smac VH3                 | SA12201L  | 1L         | Pre-packed columns feature standard connectors and are compatible with various low-to-medium pressure chromatography systems, such as ÄKTA, offering convenience for user operation. |
|                                 | SA12210L  | 10L        |  |
|                                 | SA122C11  | 1×1ml      |  |
|                                 | SA122C51  | 5×1ml      |  |
| HiPur Smac VH3                  | SA122C15  | 1×5ml      |  |
|                                 | SA122C55  | 5×5ml      |  |
| HiPur Smac VH3                  | SA122C20  | 20ml,1pcs  |  |
| HiSelect Smac VH3               | SA122C47  | 4.7ml,1pcs |  |

| Products                | Cat.No.  | Size       | Applications  |
|-------------------------|----------|------------|---|
| rProtein A-ME Beads 4FF | SA073005 | 5ml        | Primarily designed for purification of monoclonal antibodies, polyclonal antibodies, or Fc-fusion proteins, this system enables antibody elution under elevated pH conditions (pH 5.0), effectively preserving biological activity. |
|                         | SA073025 | 25ml       |   |
|                         | SA073100 | 100ml      |   |
|                         | SA073500 | 500ml      |   |
|                         | SA07301L | 1L         |   |
| HiSelect A-ME           | SA073C47 | 4.7ml,1pcs | Designed for the separation and purification of monoclonal antibodies, this product is well-suited for large-scale purification of industrial or clinically applicable antibodies.  |
| HiPur A-ME              | SA073C20 | 20ml,1pcs  |   |
| Smac SupA-ME            | SA137005 | 5ml        |   |
|                         | SA137025 | 25ml       |   |
|                         | SA137100 | 100ml      |   |
|                         | SA137500 | 500ml      |   |
|                         | SA13701L | 1L         |   |
| PreCap Smac SupA-ME     | SA137C11 | 1x1ml      |   |
|                         | SA137C51 | 5x1ml      |   |
|                         | SA137C15 | 1x5ml      |   |
| HiSelect Smac SupA-ME   | SA137C55 | 5x5ml      |   |
| HiPur Smac SupA-ME      | SA137C47 | 4.7ml,1pcs |   |
|                         | SA137C20 | 20ml,1pcs  |   |

### Monoclonal antibody purification resin

Monoclonal antibody purification products mainly refer to Protein At Beads 4FF, Protein At Beads LXPro and Protein SupAt Beads, which can purify monoclonal antibodies at relatively high flow rates. These products are designed for cleaning-in-place using 0.1-1 M NaOH, making them suitable for industrial customers aiming to efficiently produce and scale up antibody production. Key features include:

#### Product Advantages

- High dynamic binding capacity and adsorption efficiency.
- Resistance to cleaning-in-place with 0.1-1 M NaOH, ensuring a long shelf life. After 200 cleaning cycles with 0.1 M NaOH, the medium's loading capacity remains nearly unchanged. After 100 cycles with 0.5 M NaOH, the loading capacity still exceeds 80% of the initial value.
- A relatively stable directional coupling method is employed, resulting in low ligand leaching (less than 20 ng/mg IgG).
- Effective removal of host cell proteins (HCP) and DNA, achieving HCP residues below 1 µg/mg and HCDNA residues below 1 ppm, comparable to leading commercial products.
- Comprehensive documentation and data are provided to support drug submissions.

#### Product Information

| Products                            | Cat.No.   | Size  | Applications  |
|-------------------------------------|-----------|-------|---|
| Protein At Beads 4FF                | SA023005  | 5ml   | Protein At Beads 4FF is easy to pack and use. Its high stability and high flow properties make it excellent for scaling-up. |
|                                     | SA023025  | 25ml  |   |
|                                     | SA023100  | 100ml |   |
|                                     | SA023500  | 500ml |   |
|                                     | SA02301L  | 1L    |   |
|                                     | SA02310L  | 10L   |   |
| Protein At Beads 4FF Gravity Column | SA023GC01 | 1ml   |   |
|                                     | SA023GC05 | 5ml   |   |

| Products               | Cat.No.  | Size        | Applications   |
|------------------------|----------|-------------|--|
| MabCap At 4FF          | SA023C11 | 1x1ml       | The columns have the standard interface that can be adapted to the liquid chromatography systems, such as Smart GO™ and Smart PUR™.  |
|                        | SA023C51 | 5x1ml       |  |
|                        | SA023C15 | 1x5ml       |  |
|                        | SA023C55 | 5x5ml       |  |
|                        | SA023CS  | 3x1ml+1x5ml |  |
| HiSelect At 4FF        | SA023C47 | 4.7ml,1pcs  | Protein At Beads LX has a high dynamic binding capacity (DBC) and is specially designed for processing feed from high-expression cell cultures with increased antibody titers. |
| HiPur At 4FF           | SA023C20 | 20ml,1pcs   |  |
| Protein At Beads LX    | SA085005 | 5ml         |  |
|                        | SA085025 | 25ml        |  |
|                        | SA085100 | 100ml       |  |
|                        | SA085500 | 500ml       |  |
|                        | SA08501L | 1L          |  |
| MabCap At LX           | SA08510L | 10L         |  |
|                        | SA085C11 | 1x1ml       |  |
|                        | SA085C51 | 5x1ml       |  |
|                        | SA085C15 | 1x5ml       |  |
|                        | SA085C55 | 5x5ml       |  |
| HiSelect At LX         | SA085CS  | 3x1ml+1x5ml | The columns have the standard interface that can be adapted to the liquid chromatography systems, such as Smart GO™ and Smart PUR™.  |
| HiPur At LX            | SA085C47 | 4.7ml,1pcs  |  |
| Protein At Beads LXPro | SA085C20 | 20ml,1pcs   |  |
|                        | SA124005 | 5ml         |  |
|                        | SA124025 | 25ml        |  |
|                        | SA124100 | 100ml       |  |
|                        | SA124500 | 500ml       |  |
| HiSelect At LXPro      | SA12401L | 1L          |  |
|                        | SA12410L | 10L         |  |
|                        | SA124C47 | 4.7ml,1pcs  |  |
| HiPur At LXPro         | SA124C20 | 20ml,1pcs   |  |
|                        | SA124C11 | 1x1ml       | The columns have the standard interface that can be adapted to the liquid chromatography systems, such as Smart GO™ and Smart PUR™.  |
| PreCap At LXPro        | SA124C51 | 5x1ml       |  |
|                        | SA124C15 | 1x5ml       |  |
|                        | SA124C55 | 5x5ml       |  |
| Protein SupAt Beads    | SA095005 | 5ml         |  |
|                        | SA095025 | 25ml        |  |
|                        | SA095100 | 100ml       |  |
|                        | SA095500 | 500ml       |  |
| MabCap SupAt Beads     | SA09501L | 1L          |  |
|                        | SA095C11 | 1x1ml       |  |
|                        | SA095C51 | 5x1ml       |  |
|                        | SA095C15 | 1x5ml       |  |
|                        | SA095C55 | 5x5ml       |  |
| HiSelect SupAt         | SA095CS  | 3x1ml+1x5ml |  |
| HiPur SupAt            | SA095C47 | 4.7ml,1pcs  |  |
|                        | SA095C20 | 20ml,1pcs   |  |

## Protein G related products

### Universal antibody purification resin

rProtein G Beads and rProtein G Beads 4FF are utilized for the purification of antibodies. Protein G is a cell wall protein isolated from *G. streptococci*, known for its ability to bind mammalian IgG through the Fc fragment. The recombinant protein G has high-affinity binding capacity and low nonspecificity. Notably, Protein G and Protein A exhibit distinct binding properties to IgG. Compared to Protein A, Protein G demonstrates a stronger binding affinity for polyclonal antibodies from sources such as cattle, sheep, and horses. It is also capable of binding rat IgG, human IgG3, and mouse IgG1, which do not bind well with Protein A. The AbPur rProtein G kit is designed for the one-step purification of antibodies. This kit includes a prepacked column and ready-to-use buffer.

#### Example: AbCap G 4FF




Fig. 1-14 Purification of antibody by AbCap G 4FF.

### Product Information

| Products                            | Cat.No.   | Size        | Applications   |
|-------------------------------------|-----------|-------------|--|
| rProtein G Beads                    | SA016005  | 5ml         | The kit can be used to purify both monoclonal and polyclonal antibodies, and the purity can reach 95% in one step it is recommended for use at the laboratory scale. |
|                                     | SA016025  | 25ml        |  |
|                                     | SA016100  | 100ml       |  |
|                                     | SA016500  | 500ml       |  |
|                                     | SA01601L  | 1L          |  |
|                                     | SA01610L  | 10L         |  |
| rProtein G Beads Gravity Column     | SA016GC01 | 1ml         | The one-step antibody purification kit is characterized by its simplicity of operation and high purification efficiency.   |
| AbPur rProtein G Kit AbPur          | SA016K03  | 3T          |  |
| rProtein G Buffer Kit               | SA016KB03 | 3T          | rProtein G Beads 4FF can be used in FPLC and is easy to pack and use. Its high flow properties make it excellent for scaling-up.                                     |
| rProtein G Beads 4FF                | SA020005  | 5ml         |  |
|                                     | SA020025  | 25ml        |  |
|                                     | SA020100  | 100ml       |  |
|                                     | SA020500  | 500ml       |  |
|                                     | SA02001L  | 1L          |  |
| SA02010L                            | 10L       |             |  |
| rProtein G Beads 4FF Gravity Column | SA020GC01 | 1ml         | The columns have the standard interface that can be adapted to the liquid chromatography systems, such as Smart GO™ and Smart PUR™.                                  |
|                                     | SA020GC05 | 5ml         |  |
| AbCap G 4FF                         | SA020C11  | 1×1ml       |  |
|                                     | SA020C51  | 5×1ml       |  |
|                                     | SA020C15  | 1×5ml       |  |
|                                     | SA020C55  | 5×5ml       |  |
|                                     | SA020CS   | 3×1ml+1×5ml |  |
| HiSelect G 4FF HiPur G 4FF          | SA020C47  | 4.7ml,1pcs  |  |
|                                     | SA020C20  | 20ml,1pcs   |  |

## Protein A/G related products

### Universal antibody purification resin

rProtein A/G Beads 4FF is an affinity medium, designed for the purification, immunoprecipitation and Co-immunoprecipitation. The rProtein A/G is coupled to the highly cross-linked 4% agarose beads by a stable thioether linkage. The medium has higher antibody binding strength and lower protein nonspecific adsorption.

### Product Information

| Products               | Cat.No.  | Size        | Applications  |
|------------------------|----------|-------------|---|
| rProtein A/G Beads 4FF | SA032005 | 5ml         | It is mainly used for detection of antigens and antibodies such as IP and Co-IP.  |
|                        | SA032025 | 25ml        |   |
|                        | SA032100 | 100ml       |   |
|                        | SA032500 | 500ml       |   |
|                        | SA03201L | 1L          |   |
|                        | SA03210L | 10L         |   |
| AbCap A/G 4FF          | SA032C11 | 1×1ml       | The columns have the standard interface that can be adapted to the liquid chromatography systems, such as Smart GO™ and Smart PUR™. |
|                        | SA032C51 | 5×1ml       |   |
|                        | SA032C15 | 1×5ml       |   |
|                        | SA032C55 | 5×5ml       |   |
|                        | SA032CS  | 3×1ml+1×5ml |   |
| HiSelect A/G 4FF       | SA032C47 | 4.7ml,1pcs  |   |
| HiPur A/G 4FF          | SA032C20 | 20ml,1pcs   |   |

## Protein L related products

### Universal antibody purification resin

rProtein L Beads, rProtein L Beads 4FF and Smac L are universal affinity chromatography medium designed for the simple, one-step purification of various classes, subclasses, and fragments of immunoglobulins from biological samples and cell culture supernatant. Unlike Protein A and Protein G, which primarily bind to the Fc regions of immunoglobulins, Protein L targets the variable region of kappa light chains, allowing it to bind a broader range of immunoglobulin classes and subclasses. Consequently, Protein L is particularly well-suited for capturing different-sized antibody fragments, including Fab fragments, single-chain variable fragments (scFvs), and domain antibodies (dAbs).

### Product Information

| Products         | Cat.No.  | Size  | Applications  |
|------------------|----------|-------|---|
| rProtein L Beads | SA045005 | 5ml   | rProtein L is obtained by gene recombination technology. rProtein L Beads 4FF can bind to the kappa light chains of antibodies. |
|                  | SA045025 | 25ml  |   |
|                  | SA045100 | 100ml |   |
|                  | SA045500 | 500ml |   |
|                  | SA04501L | 1L    |   |
|                  | SA04510L | 10L   |   |

| Products             | Cat.No.  | Size        | Applications  |
|----------------------|----------|-------------|---|
| rProtein L Beads 4FF | SA033005 | 5ml         | rProtein L is obtained by gene recombination technology. rProtein L Beads 4FF can bind to the kappa light chains of antibodies.   |
|                      | SA033025 | 25ml        |   |
|                      | SA033100 | 100ml       |   |
|                      | SA033500 | 500ml       |   |
|                      | SA03301L | 1L          |   |
|                      | SA03310L | 10L         |   |
| AbCap L 4FF          | SA033C11 | 1×1ml       |   |
|                      | SA033C51 | 5×1ml       |   |
|                      | SA033C15 | 1×5ml       |   |
|                      | SA033C55 | 5×5ml       |   |
|                      | SA033CS  | 3×1ml+1×5ml |   |
|                      | SA033C47 | 4.7ml,1pcs  |   |
| HiSelect L 4FF       | SA033C20 | 20ml,1pcs   |   |
| Smac L               | SA100005 | 5ml         | Smac L is an affinity chromatography medium designed for the purification of both monoclonal and bispecific antibodies. It is composed of highly rigid agarose microspheres, allowing for purification at relatively high flow rates. |
|                      | SA100025 | 25ml        |   |
|                      | SA100100 | 100ml       |   |
|                      | SA100500 | 500ml       |   |
|                      | SA10001L | 1L          |   |
| PreCap Smac L        | SA100C11 | 1×1ml       |   |
|                      | SA100C51 | 5×1ml       |   |
|                      | SA100C15 | 1×5ml       |   |
|                      | SA100C55 | 5×5ml       |   |
| HiSelect Smac L      | SA100C47 | 4.7ml,1pcs  |   |
| HiPur Smac L         | SA100C20 | 20ml,1pcs   |   |

## Pre-activated resin

### Polyclonal antibody purification resin

PabPur SulfoLink Beads 4FF allow for efficient covalent immobilization of sulfhydryl-containing peptides, proteins, and other ligands for use in affinity purification. Specifically, cysteine residues at the N- or C-terminus of antigens can be effectively coupled to PabPur SulfoLink Beads, facilitating antibody purification. NHS-Activated Beads 4FF are pre-activated agarose beads to which ligands, such as proteins and peptides, can be coupled via primary amino groups. This pre-activated medium enables the preparation of affinity adsorbents capable of isolating specific substances from complex mixtures, often achieving high purity in one step.




Fig. 1-15 The process of purification of polyclonal antibody with pre-activated Beads.

## Product Information

| Products                     | Cat.No.  | Size  | Applications  |  |
|------------------------------|----------|-------|---|--|
| PabPur SulfoLink Beads       | SA018005 | 5ml   | Immobilization of sulfhydryl-containing peptides, proteins, and other ligands onto a medium for use in affinity purification.   |  |
|                              | SA018025 | 25ml  |   |  |
|                              | SA018100 | 100ml |   |  |
|                              | SA018250 | 250ml |   |  |
| PabPur SulfoLink Beads 4FF   | SA069005 | 5ml   |   |  |
|                              | SA069025 | 25ml  |   |  |
|                              | SA069100 | 100ml |   |  |
|                              | SA069250 | 250ml |   |  |
| NHS-Activated Beads 4FF      | SA039005 | 5ml   |   | NHS-Activated Beads 4FF is used for covalent immobilization of substances by the primary amino groups.   |
|                              | SA039025 | 25ml  |   |  |
|                              | SA039100 | 100ml |   |  |
|                              | SA039500 | 500ml |   |  |
| NHS-Activated LA Beads 4FF   | SA03901L | 1L    | It reduces steric hindrance between ligands, offers strong pressure resistance, and maintains stability after protein coupling, making it suitable for large-scale industrial purification.           |  |
|                              | SA083005 | 5ml   |   |  |
|                              | SA083025 | 25ml  |   |  |
|                              | SA083100 | 100ml |   |  |
|                              | SA083500 | 500ml |   |  |
| Epoxy-Activated Beads 6B     | SA08301L | 1L    | This epoxy-activated agarose microsphere is designed for direct coupling of proteins and samples containing amino, thiol, or hydroxyl groups.   |  |
|                              | SA071005 | 5ml   |   |  |
|                              | SA071025 | 25ml  |   |  |
|                              | SA071100 | 100ml |   |  |
| Epoxy-Activated Beads 4FF    | SA071500 | 500ml |   |  |
|                              | SA07101L | 1L    |   |  |
|                              | SA040005 | 5ml   |   |  |
|                              | SA040025 | 25ml  |   |  |
| Epoxy-Activated Beads 4FF    | SA040100 | 100ml |   | This resin features a free carboxyl group at the end of a 10-atom hydrophilic spacer arm, enabling direct coupling of amino-containing proteins and samples. |
|                              | SA040500 | 500ml |   |  |
|                              | SA04001L | 1L    |   |  |
|                              | SA063005 | 5ml   |   |  |
| ECH-Activated Beads 4FF      | SA063025 | 25ml  | It can be covalently coupled with carboxyl-containing molecules to prepare specialized affinity media, enabling rapid and efficient one-step purification of target substances from complex mixtures. |  |
|                              | SA063100 | 100ml |   |  |
|                              | SA063500 | 500ml |   |  |
|                              | SA06301L | 1L    |   |  |
| Amino-Activated Beads 4FF    | SA043005 | 5ml   | Suitable for direct coupling of samples containing amino or thiol groups, including polysaccharides, proteins, and peptides.  |  |
|                              | SA043025 | 25ml  |   |  |
|                              | SA043100 | 100ml |   |  |
|                              | SA043500 | 500ml |   |  |
| Aldehyde-Activated Beads 4FF | SA04301L | 1L    |   |  |
|                              | SA079005 | 5ml   |   |  |
|                              | SA079025 | 25ml  |   |  |
|                              | SA079100 | 100ml |   |  |
|                              | SA079500 | 500ml |   |  |
|                              | SA07901L | 1L    |   |  |

# Other Chromatography

## Specific groups purification resin

Blue Beads 6FF is an affinity chromatography designed for the purification of albumin, interferon, lipoproteins, and blood coagulation factors. It also binds various enzymes, including kinases, dehydrogenases, and those requiring adenylyl-containing cofactors. Based on highly cross-linked 6% agarose beads, Blue Beads 6FF has excellent chemical and mechanical stability, making it ideal for large-scale purification of albumin and interferon.




Fig.1-16 One step purification of HSA from human serum.

Benzamidine Beads 4FF is based on highly cross-linked 4% agarose beads with p-aminobenzamidine covalently coupled as the ligand. As a broad-spectrum inhibitor of serine proteases, Benzamidine Beads 4FF are ideal for purifying serine proteinases, including trypsin, trypsin-like proteases, urokinase, and prekallikrein.




Fig.1-17 The affinity trypsin of chymotrypsin by Benzamidine Beads 4FF.

## Product Information

| Products                         | Cat.No.   | Size       | Applications   |
|----------------------------------|-----------|------------|--|
| Heparin Beads 6FF                | SA024005  | 5ml        | Heparin Beads 6FF is designed for the purification of proteins such as antithrombin III, coagulation factors, and other plasma proteins, as well as DNA-binding proteins, lipoproteins, protein synthesis factors, nucleic acid-acting enzymes, and steroid receptors. |
|                                  | SA024025  | 25ml       |  |
|                                  | SA024100  | 100ml      |  |
|                                  | SA024500  | 500ml      |  |
|                                  | SA02401L  | 1L         |  |
| Heparin Beads 6FF Gravity Column | SA024GC01 | 1ml        |  |
|                                  | SA024GC05 | 5ml        |  |
| PreCap Heparin                   | SA024C11  | 1×1ml      | The columns have the standard interface that can be adapted to the liquid chromatography systems, such as Smart GO™ and Smart PUR™.  |
|                                  | SA024C51  | 5×1ml      |  |
|                                  | SA024C15  | 1×5ml      |  |
|                                  | SA024C55  | 5×5ml      |  |
| HiSelect Heparin 6FF             | SA024C47  | 4.7ml,1pcs |  |
| HiPur Heparin 6FF                | SA024C20  | 20ml,1pcs  |  |
| Blue Beads 6FF                   | SA027005  | 5ml        | It can purify albumin, interferon, lipoproteins and blood coagulation factors and most enzymes requiring adenylyl-containing cofactors. It is an ideal medium for purification of albumin and interferon at process scale.   |
|                                  | SA027025  | 25ml       |  |
|                                  | SA027100  | 100ml      |  |
|                                  | SA027500  | 500ml      |  |
|                                  | SA02701L  | 1L         |  |
| PreCap Blue                      | SA027C11  | 1×1ml      | The columns have the standard interface that can be adapted to the liquid chromatography systems, such as Smart GO™ and Smart PUR™.  |
|                                  | SA027C51  | 5×1ml      |  |
|                                  | SA027C15  | 1×5ml      |  |
|                                  | SA027C55  | 5×5ml      |  |
| HiSelect Blue 6FF                | SA027C47  | 4.7ml,1pcs |  |
| HiPur Blue 6FF                   | SA027C20  | 20ml,1pcs  |  |
| Smac Blue                        | SA128005  | 5ml        | Designed for the purification of albumin, interferons, nucleotide-requiring enzymes, α2-macroglobulin, coagulation factors, and other related biomolecules.  |
|                                  | SA128025  | 25ml       |  |
|                                  | SA128100  | 100ml      |  |
|                                  | SA128500  | 500ml      |  |
|                                  | SA12801L  | 1L         |  |
|                                  | SA12810L  | 10L        |  |
| PreCap Smac Blue                 | SA128C11  | 1×1ml      | The columns have the standard interface that can be adapted to the liquid chromatography systems, such as Smart GO™ and Smart PUR™.  |
|                                  | SA128C51  | 5×1ml      |  |
|                                  | SA128C15  | 1×5ml      |  |
| HiSelect Smac Blue               | SA128C47  | 4.7ml,1pcs |  |
| HiPur Smac Blue                  | SA128C20  | 20ml,1pcs  |  |
| Con A Beads 4FF                  | SA028005  | 5ml        | Con A Beads 4FF can detect and purify glycoproteins, membrane proteins, glycolipids, polysaccharides, membrane vesicles with mannoside or glucoside residues, IgM, hormones and lipoproteins.  |
|                                  | SA028025  | 25ml       |  |
|                                  | SA028100  | 100ml      |  |
|                                  | SA028500  | 500ml      |  |
|                                  | SA02801L  | 1L         |  |
| Con A Beads 4FF Gravity Column   | SA028GC01 | 1ml        |  |
|                                  | SA028GC05 | 5ml        |  |

| Products                              | Cat.No.   | Size        | Applications  |
|---------------------------------------|-----------|-------------|---|
| PreCap Con A                          | SA028C11  | 1x1ml       | The columns have the standard interface that can be adapted to the liquid chromatography systems, such as Smart GO™ and Smart PUR™.   |
|                                       | SA028C51  | 5x1ml       |   |
|                                       | SA028C15  | 1x5ml       |   |
|                                       | SA028C55  | 5x5ml       |   |
|                                       | SA028CS   | 3x1ml+1x5ml |   |
| HiSelect Con A 4FF                    | SA028C47  | 4.7ml,1pcs  |   |
| HiPur Con A 4FF                       | SA028C20  | 20ml,1pcs   |   |
| Endotoxin Removal Beads               | SA031005  | 5ml         | Polymyxin B is coupled to 4% agarose beads, enabling effective endotoxin removal from biological samples. The ligand of this product is PMB affinity ligand.                                |
|                                       | SA031025  | 25ml        |   |
|                                       | SA031100  | 100ml       |   |
|                                       | SA031500  | 500ml       |   |
|                                       | SA03101L  | 1L          |   |
| Endotoxin Removal Bead Gravity Column | SA031GC01 | 1ml         |   |
|                                       | SA031GC05 | 5ml         |   |
| Endotoxin Removal Kit                 | SA031K03  | 3T          | It contains gravity columns and buffers for endotoxin removal.  |
| LCA Beads 4FF                         | SA097005  | 5ml         | LCA Beads 4FF can detect and purify glycoproteins, membrane proteins, glycolipids, polysaccharides, membrane vesicles with mannoside or glucoside residues, IgM, hormones and lipoproteins. |
|                                       | SA097025  | 25ml        |   |
|                                       | SA097100  | 100ml       |   |
|                                       | SA097500  | 500ml       |   |
|                                       | SA09701L  | 1L          |   |
| PreCap LCA                            | SA097C11  | 1x1ml       |   |
|                                       | SA097C51  | 5x1ml       |   |
|                                       | SA097C15  | 1x5ml       |   |
|                                       | SA097C55  | 5x5ml       |   |
| Smac LCA                              | SA097CS   | 3x1ml+1x5ml |   |
|                                       | SA127005  | 5ml         |   |
|                                       | SA127025  | 25ml        |   |
|                                       | SA127100  | 100ml       |   |
|                                       | SA127500  | 500ml       |   |
|                                       | SA12701L  | 1L          |   |
|                                       | SA12710L  | 10L         |   |
| Benzamidine Beads 4FF                 | SA044005  | 5ml         |   |
|                                       | SA044025  | 25ml        |   |
|                                       | SA044100  | 100ml       |   |
|                                       | SA044500  | 500ml       |   |
|                                       | SA04401L  | 1L          |   |
|                                       | SA04410L  | 10L         |   |
| PreCap Benzamidine 4FF                | SA044C11  | 1x1ml       | The columns have the standard interface that can be adapted to the liquid chromatography systems, such as Smart GO™ and Smart PUR™.   |
|                                       | SA044C51  | 5x1ml       |   |
|                                       | SA044C15  | 1x5ml       |   |
|                                       | SA044C55  | 5x5ml       |   |
|                                       | SA044CS   | 3x1ml+1x5ml |   |
| HiSelect Benzamidine 4FF              | SA044C47  | 4.7ml,1pcs  |   |
| HiPur Benzamidine 4FF                 | SA044C20  | 20ml,1pcs   |   |

| Products                    | Cat.No.  | Size       | Applications   |
|-----------------------------|----------|------------|--|
| Iminobiotin Beads 6FF       | SA059005 | 5ml        | It enables one-step purification of streptavidin.  |
|                             | SA059025 | 25ml       |  |
|                             | SA059100 | 100ml      |  |
|                             | SA059500 | 500ml      |  |
|                             | SA05901L | 1L         |  |
|                             | SA05910L | 10L        |  |
| HiSelect Iminobiotin 6FF    | SA059C47 | 4.7ml,1pcs |  |
| HiPur Iminobiotin 6FF       | SA059C20 | 20ml,1pcs  |  |
| Gelatin Beads 4FF           | SA054005 | 5ml        | It enables one-step purification of fibronectin.   |
|                             | SA054025 | 25ml       |  |
|                             | SA054100 | 100ml      |  |
|                             | SA054500 | 500ml      |  |
|                             | SA05401L | 1L         |  |
| Boric Acid Beads 4FF        | SA0410L  | 10L        | It can be used to purify carbohydrates, nucleosides, nucleotides, nucleic acids and enzymes, especially for the glycoproteins.                           |
|                             | SA057005 | 5ml        |  |
|                             | SA057025 | 25ml       |  |
|                             | SA057100 | 100ml      |  |
|                             | SA057500 | 500ml      |  |
|                             | SA05701L | 1L         |  |
| Lysine Beads 4FF            | SA05710L | 10L        | It can be used to purify fibrinolytic activator, fibrinolytic progenitor and ribonucleic acid.   |
|                             | SA046005 | 5ml        |  |
|                             | SA046025 | 25ml       |  |
|                             | SA046100 | 100ml      |  |
|                             | SA046500 | 500ml      |  |
| Smac PlasmidSelect          | SA04601L | 1L         |  |
|                             | SA04610L | 10L        |  |
|                             | SA104005 | 5ml        |  |
|                             | SA104025 | 25ml       |  |
| PreCap PlasmidSelect        | SA104100 | 100ml      | This thiophilic aromatic chromatography medium selectively separates covalently closed circular supercoiled plasmid DNA from its open-circular isoforms. |
|                             | SA104500 | 500ml      |  |
|                             | SA10401L | 1L         |  |
|                             | SA104C11 | 1x1 ml     |  |
| HiSelect Smac PlasmidSelect | SA104C51 | 5x1 ml     |  |
|                             | SA104C15 | 1x5 ml     |  |
|                             | SA104C55 | 5x5 ml     |  |
| HiPur Smac PlasmidSelect    | SA104C47 | 4.7ml,1pcs |  |
| Plasmid Purification Beads  | SA104C20 | 20ml,1pcs  | It enables rapid, one-step purification of plasmid.  |
|                             | SI014030 | 30ml       |  |
|                             | SI014100 | 100ml      |  |
|                             | SI014300 | 300ml      |  |
|                             | SI014500 | 500ml      |  |
|                             | SI01401L | 1L         |  |

| Products               | Cat.No.  | Size       | Applications   |
|------------------------|----------|------------|--|
| AAV Affinity Beads 4FF | SA096005 | 5ml        | AAV Affinity Beads 4FF is an affinity chromatography medium for the purification of adeno-associated virus (AAV) with high affinity for several subtypes of AAV. |
|                        | SA096025 | 25ml       |  |
|                        | SA096100 | 100ml      |  |
|                        | SA09601L | 1L         |  |
| Smac DeVirS            | SA226005 | 5 ml       | For virus capture and intermediate purification.   |
|                        | SA226025 | 25 ml      |  |
|                        | SA226100 | 100 ml     |  |
|                        | SA226500 | 500 ml     |  |
|                        | SA22601L | 1 L        |  |
|                        | SA22610L | 10 L       |  |
|                        | SA226C11 | 1x1 ml     |  |
| PreCap Smac DeVirS     | SA226C51 | 5x1 ml     |  |
|                        | SA226C15 | 1x5 ml     |  |
|                        | SA226C55 | 5x5 ml     |  |
| HiSelect Smac DeVirS   | SA226C47 | 4.7ml,1pcs |  |
| HiPur Smac DeVirS      | SA226C20 | 20ml,1pcs  |  |

## 1.2 Ion Exchange purification resin

Smart-Lifesciences provides ion exchange resins with four ligands (DEAE, CM, SP, and Q) and two substrates. Available product formats include both resin and pre-packed columns.



Fig.1-18 The structure diagram of four ion exchange chromatography.

## Strong cation exchange chromatography resin

### Product Information

| Products            | Cat.No.  | Size       | Applications  |
|---------------------|----------|------------|---|
| SP Beads 6FF        | SI003025 | 25ml       | A strong cation exchange medium designed for separation and purification of positively charged proteins at both laboratory and industrial scales.   |
|                     | SI003100 | 100ml      |   |
|                     | SI003500 | 500ml      |   |
|                     | SI00301L | 1L         |   |
| lexCap SP 6FF       | SI00310L | 10L        | Pre-packed with SP Beads 6FF, these columns feature standard connectors and are compatible with various low-to-medium pressure chromatography systems, such as ÄKTA, offering convenience for user operation. |
|                     | SI003C11 | 1x1ml      |   |
|                     | SI003C51 | 5x1ml      |   |
|                     | SI003C15 | 1x5ml      |   |
| HiSelect SP 6FF     | SI003C55 | 5x5ml      |   |
|                     | SI003C47 | 4.7ml,1pcs |   |
| HiPur SP 6FF        | SI003C20 | 20ml,1pcs  |   |
| SP Beads 6FF        | SI011050 | 50ml       | Designed for downstream separation and purification of proteins, nucleic acids, and peptides in biopharmaceutical and bioengineering applications.  |
|                     | SI011075 | 75ml       |   |
|                     | SI011250 | 250ml      |   |
|                     | SI01101L | 1L         |   |
| SP Big Beads        | SI0110L  | 10L        |   |
|                     | SI043005 | 5ml        |   |
|                     | SI043025 | 25ml       |   |
|                     | SI043100 | 100ml      |   |
| Smac SP             | SI043500 | 500ml      | A strong cation exchange medium suitable for purification of positively charged proteins at both laboratory and industrial scales.  |
|                     | SI04301L | 1L         |   |
|                     | SI04310L | 10L        |   |
|                     | SI027025 | 25ml       |   |
| lexCap Smac SP      | SI027100 | 100ml      | Pre-packed with Smac SP, these columns feature standard connectors and are compatible with various low-to-medium pressure chromatography systems, such as ÄKTA, for ease of use.                              |
|                     | SI027500 | 500ml      |   |
|                     | SI02701L | 1L         |   |
|                     | SI02710L | 10L        |   |
| HiSelect Smac SP    | SI027C11 | 1x1ml      |   |
|                     | SI027C51 | 5x1ml      |   |
| HiPur Smac SP       | SI027C15 | 1x5ml      |   |
|                     | SI027C55 | 5x5ml      |   |
| Smac SP 40          | SI027C47 | 4.7ml,1pcs | A strong cation exchange medium for separation and purification of positively charged proteins at laboratory and industrial scales.   |
|                     | SI027C20 | 20ml,1pcs  |   |
|                     | SI028025 | 25ml       |   |
|                     | SI028100 | 100ml      |   |
| lexCap Smac SP 40   | SI028500 | 500ml      | Pre-packed with Smac SP 40, these columns feature standard connectors and are compatible with various low-to-medium pressure chromatography systems, such as ÄKTA, ensuring operational convenience.          |
|                     | SI02801L | 1L         |   |
|                     | SI02810L | 10L        |   |
|                     | SI028C11 | 1x1ml      |   |
| HiSelect Smac SP 40 | SI028C51 | 5x1ml      |   |
|                     | SI028C15 | 1x5ml      |   |
| HiPur Smac SP 40    | SI028C55 | 5x5ml      |   |
|                     | SI028C47 | 4.7ml,1pcs |   |
| Smac SP 40 Plus     | SI028C20 | 20ml,1pcs  | A strong cation exchange medium characterized by high capacity and high resolution, suitable for high-throughput intermediate purification of various biomolecules.   |
|                     | SI044025 | 25ml       |   |
|                     | SI044100 | 100ml      |   |
|                     | SI044500 | 500ml      |   |
|                     | SI04401L | 1L         |   |
|                     | SI04410L | 10L        |   |

## Weak cation exchange chromatography resin

### Product Information

| Products        | Cat.No.  | Size       | Applications  |
|-----------------|----------|------------|---|
| CM Beads 6FF    | SI007025 | 25ml       | A weak cation exchange medium designed for separation and purification of positively charged proteins at both laboratory and industrial scales.   |
|                 | SI007100 | 100ml      |   |
|                 | SI007500 | 500ml      |   |
|                 | SI00701L | 1L         |   |
| lexCap CM 6FF   | SI00710L | 10L        | Pre-packed with CM Beads 6FF, these columns feature standard connectors and are compatible with various low-to-medium pressure chromatography systems, such as ÄKTA, offering convenience for user operation. |
|                 | SI007C11 | 1×1ml      |   |
|                 | SI007C51 | 5×1ml      |   |
|                 | SI007C15 | 1×5ml      |   |
| HiSelect CM 6FF | SI007C55 | 5×5ml      |   |
|                 | SI007C47 | 4.7ml,1pcs |   |
| HiPur CM 6FF    | SI007C20 | 20ml,1pcs  |   |

## Strong anion exchange chromatography resin

### Product Information

| Products           | Cat.No.  | Size       | Applications   |
|--------------------|----------|------------|--|
| Q Beads 6FF        | SI001025 | 25ml       | A strong anion exchange medium designed for separation and purification of negatively charged proteins at both laboratory and industrial scales.   |
|                    | SI001100 | 100ml      |  |
|                    | SI001500 | 500ml      |  |
|                    | SI00101L | 1L         |  |
| lexCap Q 6FF       | SI00110L | 10L        | Pre-packed with QBeads 6FF, these columns feature standard connectors and are compatible with various low-to-medium pressure chromatography systems, such as ÄKTA, for ease of use.                        |
|                    | SI001C11 | 4×1ml      |  |
|                    | SI001C51 | 5×1ml      |  |
|                    | SI001C15 | 1×5ml      |  |
| HiSelect Q 6FF     | SI001C55 | 5×5ml      |  |
| HiPur Q 6FF        | SI001C47 | 4.7ml,1pcs |  |
| Q Beads HP         | SI001C20 | 20ml,1pcs  | A strong anion exchange medium suitable for separation and purification of negatively charged proteins at both laboratory and industrial scales.   |
|                    | SI010050 | 50ml       |  |
|                    | SI010075 | 75ml       |  |
|                    | SI010250 | 250ml      |  |
| Smac Q             | SI01001L | 1L         |  |
|                    | SI01010L | 10L        |  |
|                    | SI019025 | 25ml       |  |
|                    | SI019100 | 100ml      |  |
| lexCap Smac Q      | SI019500 | 500ml      | Pre-packed with Smac Q, these columns feature standard connectors and are compatible with various low-to-medium pressure chromatography systems, such as ÄKTA, ensuring operational convenience.           |
|                    | SI01901L | 1L         |  |
|                    | SI01910L | 10L        |  |
|                    | SI019C11 | 1×1ml      |  |
| HiSelect Smac Q    | SI019C51 | 5×1ml      |  |
|                    | SI019C15 | 1×5ml      |  |
| HiPur Smac Q       | SI019C55 | 5×5ml      |  |
| Smac Q 40          | SI019C47 | 4.7ml,1pcs | A strong anion exchange medium designed for separation and purification of positively charged proteins at both laboratory and industrial scales.   |
|                    | SI035025 | 25ml       |  |
|                    | SI035100 | 100ml      |  |
|                    | SI035500 | 500ml      |  |
| lexCap Smac Q 40   | SI03501L | 1L         | Pre-packed with Smac Q 40, these columns feature standard connectors and are compatible with various low-to-medium pressure chromatography systems, such as ÄKTA, offering convenience for user operation. |
|                    | SI03510L | 10L        |  |
|                    | SI035C11 | 1×1ml      |  |
|                    | SI035C51 | 5×1ml      |  |
| HiSelect Smac Q 40 | SI035C15 | 1×5ml      |  |
|                    | SI035C55 | 5×5ml      |  |
| HiPur Smac Q 40    | SI035C47 | 4.7ml,1pcs |  |
|                    | SI035C20 | 20ml,1pcs  |  |

## Weak anion exchange chromatography resin

### Product Information

| Products           | Cat.No.  | Size       | Applications  |
|--------------------|----------|------------|---|
| DEAE Beads 6FF     | SI005025 | 25ml       | A weak anion exchange medium designed for separation and purification of negatively charged proteins at both laboratory and industrial scales.  |
|                    | SI005100 | 100ml      |   |
|                    | SI005500 | 500ml      |   |
|                    | SI00501L | 1L         |   |
| lexCap DEAE 6FF    | SI00510L | 10L        | Pre-packed with DEAE Beads 6FF, these columns feature standard connectors and are compatible with various low-to-medium pressure chromatography systems, such as ÄKTA, offering convenience for user operation. |
|                    | SI005C11 | 4×1ml      |   |
|                    | SI005C51 | 5×1ml      |   |
|                    | SI005C15 | 1×5ml      |   |
| HiSelect DEAE 6FF  | SI005C55 | 5×5ml      |   |
| HiPur DEAE 6FF     | SI007C47 | 4.7ml,1pcs |   |
| DEAE Smartdex A-25 | SI007C20 | 20ml,1pcs  | A dextran-based weak anion exchange medium. Widely used in downstream separation and purification of proteins, nucleic acids, and peptides in biopharmaceutical and bioengineering applications.                |
|                    | SI0091   | 25g        |   |
|                    | SI0092   | 100g       |   |
|                    | SI0093   | 500g       |   |
|                    | SI0094   | 1Kg        |   |
|                    | SI0095   | 10Kg       |   |
|                    | SI0096   | 25Kg       |   |

## 1.3 Hydrophobic Interaction purification resin

Hydrophobic interaction chromatography consists of four ligands (butyl, octyl, phenyl and butyl-S) immobilized on highly cross-linked agarose beads. It can be used in multiple formats, including conventional gravity-flow chromatography and purification system. PreCap columns are pre-packed columns containing this medium.

### Butyl

### Product Information

| Products           | Cat.No.  | Size       | Applications   |
|--------------------|----------|------------|--|
| Butyl Beads 4FF    | SH001025 | 25ml       | Butyl Beads 4FF is an aliphatic hydrophobic interaction medium and has the weaker hydrophobicity to the substances.  |
|                    | SH001100 | 100ml      |  |
|                    | SH001500 | 500ml      |  |
|                    | SH00101L | 1L         |  |
| PreCap Butyl 4FF   | SH00101L | 1L         | The columns have the standard interface that can be adapted to the liquid chromatography systems, such as Smart GO™ and Smart PUR™.  |
|                    | SH001C11 | 1×1ml      |  |
|                    | SH001C51 | 5×1ml      |  |
|                    | SH001C15 | 1×5ml      |  |
| HiSelect Butyl 4FF | SH001C55 | 5×5ml      |  |
| HiPur Butyl 4FF    | SH001C47 | 4.7ml,1pcs |  |
| Butyl Beads HF     | SH001C20 | 20ml,1pcs  | Suitable for the separation of biomolecules with aliphatic ligands, the highly rigid agarose microspheres can withstand higher flow rates and provide enhanced chemical stability. |
|                    | SH016025 | 25ml       |  |
|                    | SH016100 | 100ml      |  |
|                    | SH016500 | 500ml      |  |
|                    | SH01601L | 1L         |  |
|                    | SH01610L | 10L        |  |

| Products               | Cat.No.  | Size       | Applications   |   |
|------------------------|----------|------------|--|---|
| HiPur Butyl HF         | SH016C20 | 20ml,1pcs  | Packed with Butyl Beads HF, the prepacked column features standard fittings for compatibility with a wide range of low-to-medium pressure chromatography systems, simplifying customer operation.                                  |   |
|                        | SH016C11 | 1x1ml      |  |   |
|                        | SH016C51 | 5x1ml      |  |   |
|                        | SH016C15 | 1x5ml      |  |   |
| PreCap Butyl HF        | SH016C55 | 5x5ml      | Engineered for the separation of biomolecules with aliphatic ligands, this chromatography resin delivers high binding capacity and small particle size, ensuring reliable performance for industrial-scale purification processes. |   |
|                        | SH016C47 | 4.7ml,1pcs |  |   |
| HiSelect Butyl HF      | SH014025 | 25ml       | Packed with Smac Butyl 40, the column features standard fittings for seamless compatibility with diverse low-to-medium pressure chromatography systems, streamlining user operation.   |   |
|                        | SH014100 | 100ml      |  |   |
|                        | SH014500 | 500ml      |  |   |
|                        | SH01401L | 1L         |  |   |
|                        | SH01410L | 10L        |  |   |
| Smac Butyl 40          | SH014C47 | 4.7ml,1pcs | Designed for binding and elution of relatively hydrophobic molecules under moderately low salt conditions, suitable for separation or purification of various biomolecules at laboratory or process scales.                        |   |
| HiSelect Smac Butyl 40 | SH014C20 | 20ml,1pcs  |  |   |
|                        | SH014C11 | 1x1ml      |  |   |
| HiPur Smac Butyl 40    | SH014C51 | 5x1ml      |  |   |
|                        | SH014C15 | 1x5ml      |  |   |
|                        | SH014C55 | 5x5ml      |  |   |
| PreCap Butyl 40        | SH012005 | 5ml        |  |   |
|                        | SH012025 | 25ml       |  |   |
|                        | SH012100 | 100ml      |  |   |
|                        | SH01201L | 1L         |  |   |
| Butyl-S Beads 6FF      | SH01210L | 10L        |  | The columns have the standard interface that can be adapted to the liquid chromatography systems, such as Smart GO™ and Smart PUR™. |
|                        | SH012C11 | 1x1m       |  |   |
|                        | SH012C51 | 5x1ml      |  |   |
|                        | SH012C15 | 1x5ml      |  |   |
| PreCap Butyl-S 6FF     | SH012C55 | 5x5ml      | Suitable for the separation and purification of proteins with aromatic ligands, the highly rigid agarose microspheres withstand higher flow rates and provide enhanced chemical stability.   |   |
|                        | SH012C47 | 4.7ml,1pcs |  |   |
|                        | SH012C20 | 20ml,1pcs  |  |   |
| HiSelect Butyl-S       |          |            | Packed with Phenyl Beads HF (High Sub), the prepacked column features standard fittings for compatibility with a wide range of low-to-medium pressure chromatography systems, enhancing user-friendliness.                         |   |
| HiPur Butyl-S          |          |            |  |   |

## Octyl

### Product Information

| Products           | Cat.No.  | Size       | Applications  |
|--------------------|----------|------------|---|
| Octyl Beads 4FF    | SH003025 | 25ml       | Octyl Beads 4FF is an aliphatic hydrophobic interaction medium and has the medium hydrophobicity to the substances.   |
|                    | SH003100 | 100ml      |   |
|                    | SH003500 | 500ml      |   |
|                    | SH00301L | 1L         |   |
| PreCap Octyl 4FF   | SH003C11 | 1x1ml      | The columns have the standard interface that can be adapted to the liquid chromatography systems, such as Smart GO™ and Smart PUR™.   |
|                    | SH003C51 | 5x1ml      |   |
|                    | SH003C15 | 1x5ml      |   |
|                    | SH003C55 | 5x5ml      |   |
| HiSelect Octyl 4FF | SH003C47 | 4.7ml,1pcs | Packed with Phenyl Beads HF (Low Sub), the pre-packed column features standard fittings for compatibility with diverse low-to-medium pressure chromatography systems, simplifying user operation. |
| HiPur Octyl 4FF    | SH003C20 | 20ml,1pcs  |   |

## Phenyl

### Product Information

| Products                    | Cat.No.  | Size       | Applications  |  |
|-----------------------------|----------|------------|---|--|
| Phenyl Beads 6FF (Low Sub)  | SH004025 | 25ml       | Phenyl Beads 6FF is suitable for separation of hydrophobic proteins. 1)According to the selectivity and binding capacity required for separation, there are two levels of density of the ligand phenyl. 2)The medium with low substitution (LS) and high substitution (HS) are composed of 90 um highly cross-linked agarose beads. 3)Phenyl is coupled to the uncharged agarose beads by ether bond. 4)The medium has high resolution and binding capacity. The prepacked columns are ready for use. |  |
|                             | SH004100 | 100ml      |   |  |
|                             | SH004500 | 500ml      |   |  |
|                             | SH00401L | 1L         |   |  |
| PreCap Phenyl LS 6FF        | SH004C11 | 1x1ml      |   | PreCap Select content of four 1ml prepacked columns, including PreCap Butyl 4FF, PreCap Octyl 4FF, PreCap Phenyl LS 6FF and PreCap Phenyl HS 6FF.  |
|                             | SH004C51 | 5x1ml      |   |  |
|                             | SH004C15 | 1x5ml      |   |  |
|                             | SH004C55 | 5x5ml      |   |  |
| HiSelect Phenyl 6FF LS      | SH004C47 | 4.7ml,1pcs |   | Suitable for the separation and purification of proteins with aromatic ligands, the highly rigid agarose microspheres withstand higher flow rates and provide enhanced chemical stability. |
| HiPur Phenyl 6FF LS         | SH004C20 | 20ml,1pcs  |   |  |
| Phenyl Beads 6FF (High Sub) | SH006025 | 25ml       |   |  |
|                             | SH006100 | 200ml      |   |  |
|                             | SH006500 | 500ml      |   |  |
|                             | SH00601L | 1L         |   |  |
| PreCap Phenyl HS 6FF        | SH006C11 | 1x1ml      | Packed with Phenyl Beads HF (High Sub), the prepacked column features standard fittings for compatibility with a wide range of low-to-medium pressure chromatography systems, enhancing user-friendliness.  |  |
|                             | SH006C51 | 5x1ml      |   |  |
|                             | SH006C15 | 1x5ml      |   |  |
|                             | SH006C55 | 5x5ml      |   |  |
| HiSelect Phenyl 6FF HS      | SH006C47 | 4.7ml,1pcs | Suitable for the separation and purification of proteins with aromatic ligands, the highly rigid agarose microspheres withstand higher flow rates and provide enhanced chemical stability.  |  |
| HiPur Phenyl 6FF HS         | SH006C20 | 20ml,1pcs  |   |  |
| PreCap Select               | SH015025 | 25ml       |   |  |
|                             | SH015100 | 100ml      |   |  |
|                             | SH015500 | 500ml      |   |  |
|                             | SH01501L | 1L         |   |  |
| Phenyl Beads HF(High Sub)   | SH01510L | 10L        | Packed with Phenyl Beads HF (Low Sub), the pre-packed column features standard fittings for compatibility with diverse low-to-medium pressure chromatography systems, simplifying user operation.   |  |
|                             | SH015C47 | 4.7ml,1pcs |   |  |
|                             | SH015C20 | 20ml,1pcs  |   |  |
|                             | SH015C11 | 1x1ml      |   |  |
| HiSelect Phenyl HF HS       | SH015C51 | 5x1ml      | Optimized for the separation and purification of proteins bearing aromatic ligands, this resin delivers high binding capacity and small particle size, making it a reliable solution for industrial-scale purification.   |  |
|                             | SH015C15 | 1x5ml      |   |  |
|                             | SH015C55 | 5x5ml      |   |  |
|                             | SH018025 | 25ml       |   |  |
| PreCap Phenyl HF (HS)       | SH018100 | 100ml      | Packed with smac phenyl 40, the column features standard fittings for seamless compatibility with diverse low-to-medium pressure chromatography systems, streamlining user operation.   |  |
|                             | SH018500 | 500ml      |   |  |
|                             | SH01801L | 1L         |   |  |
|                             | SH01810L | 10L        |   |  |
| Phenyl Beads HF(Low Sub)    | SH018C47 | 4.7ml,1pcs | The columns have the standard interface that can be adapted to the liquid chromatography systems, such as Smart GO™ and Smart PUR™.   |  |
|                             | SH018C20 | 20ml,1pcs  |   |  |
|                             | SH018C11 | 1x1ml      |   |  |
|                             | SH018C51 | 5x1ml      |   |  |
| HiSelect Phenyl HF LS       | SH018C15 | 1x5ml      | The columns have the standard interface that can be adapted to the liquid chromatography systems, such as Smart GO™ and Smart PUR™.   |  |
|                             | SH018C55 | 5x5ml      |   |  |
|                             | SH013025 | 25ml       |   |  |
|                             | SH013100 | 100ml      |   |  |
| HiPur Phenyl HF LS          | SH013500 | 500ml      | Optimized for the separation and purification of proteins bearing aromatic ligands, this resin delivers high binding capacity and small particle size, making it a reliable solution for industrial-scale purification.   |  |
|                             | SH01301L | 1L         |   |  |
|                             | SH01310L | 10L        |   |  |
|                             | SH013C47 | 4.7ml,1pcs |   |  |
| PreCap Phenyl HF (LS)       | SH013C20 | 20ml,1pcs  | Packed with smac phenyl 40, the column features standard fittings for seamless compatibility with diverse low-to-medium pressure chromatography systems, streamlining user operation.   |  |
|                             | SH013C11 | 1x1ml      |   |  |
|                             | SH013C51 | 5x1ml      |   |  |
|                             | SH013C15 | 1x5ml      |   |  |
| Smac Phenyl 40              | SH013C55 | 5x5ml      | The columns have the standard interface that can be adapted to the liquid chromatography systems, such as Smart GO™ and Smart PUR™.   |  |
|                             | SH013025 | 25ml       |   |  |
|                             | SH013100 | 100ml      |   |  |
|                             | SH013500 | 500ml      |   |  |
| HiSelect Smac Phenyl 40     | SH01301L | 1L         | The columns have the standard interface that can be adapted to the liquid chromatography systems, such as Smart GO™ and Smart PUR™.   |  |
|                             | SH01310L | 10L        |   |  |
|                             | SH013C47 | 4.7ml,1pcs |   |  |
|                             | SH013C20 | 20ml,1pcs  |   |  |
| HiPur Smac Phenyl 40        | SH013C11 | 1x1ml      | The columns have the standard interface that can be adapted to the liquid chromatography systems, such as Smart GO™ and Smart PUR™.   |  |
|                             | SH013C51 | 5x1ml      |   |  |
|                             | SH013C15 | 1x5ml      |   |  |
|                             | SH013C55 | 5x5ml      |   |  |
| PreCap Phenyl 40            | SH013025 | 25ml       | The columns have the standard interface that can be adapted to the liquid chromatography systems, such as Smart GO™ and Smart PUR™.   |  |
|                             | SH013100 | 100ml      |   |  |
|                             | SH013500 | 500ml      |   |  |
|                             | SH01301L | 1L         |   |  |

## 1.4 Size Exclusion Chromatography

Size exclusion chromatography (SEC), also known as gel filtration or molecular sieve chromatography, plays a vital role in separating molecules based on size as they pass through an SEC-packed column. Separation is mainly based on the size and shape of the samples, with smaller molecules being more delayed than larger ones.

### Smartarose

Smartarose is a gel filtration medium based on agarose beads, providing a broad range of separation capabilities. It consists of a porous matrix of spherical particles with excellent chemical stability and low non-specificity. By controlling the agarose concentration and degree of cross-linking, Smartarose Beads with different pore sizes can be prepared, enabling the effective separation of proteins with different molecular weights.

| Products         | Separation Range (Da) (globulin) | Recommended Flow rate | Working pH |
|------------------|----------------------------------|-----------------------|------------|
| Smartarose 4B    | 70000-20×10 <sup>6</sup>         | <30                   | 4-9        |
| Smartarose 6B    | 10000-4×10 <sup>6</sup>          | <30                   | 4-9        |
| Smartarose CL-4B | 70000-20×10 <sup>6</sup>         | <50                   | 3-13       |
| Smartarose CL-6B | 10000-4×10 <sup>6</sup>          | <80                   | 3-13       |
| Smartarose 4FF   | 40000-30×10 <sup>6</sup>         | 150-250               | 2-12       |
| Smartarose 6FF   | 10000-4×10 <sup>6</sup>          | 200-400               | 2-12       |

### Product Information

| Products                | Cat.No.   | Size       | Applications   |  |
|-------------------------|-----------|------------|--|--|
| Smartarose 4B           | SEC0060   | 25ml       | Smartarose 4B and 6B are gel filtration media based on agarose beads. They exhibit low nonspecificity and can be used for the separation of proteins, polysaccharides, and substances with different molecular weights.  |  |
|                         | SEC0061   | 100ml      |  |  |
|                         | SEC0062   | 500ml      |  |  |
|                         | SEC0063   | 1L         |  |  |
|                         | SEC0064   | 10L        |  |  |
| Smartarose 6B           | SEC0070   | 25ml       |  |  |
|                         | SEC0071   | 100 ml     |  |  |
|                         | SEC0072   | 500 ml     |  |  |
|                         | SEC0073   | 1L         |  |  |
| Smartarose CL-4B        | SEC0074   | 10L        |  |  |
|                         | SEC0080   | 25ml       | Smartarose CL-4B and CL-6B are cross-linked derivatives of Smartarose 4B and 6B, offering greater chemical and physical stability than the original products series. These CL variants can withstand certain concentrations of organic solvents, and they can be cleaned in place (CIP) with sodium hydroxide or sterilized. They are commonly used for the separation of large proteins, peptides, and polysaccharides, as well as for molecules that cannot be solubilized or agglomerated in aqueous solutions, and for molecular weight determination. |  |
|                         | SEC0081   | 100ml      |  |  |
|                         | SEC0082   | 500ml      |  |  |
| SEC0083                 | 1L        |            |  |  |
| Smartarose CL-6B        | SEC0084   | 10L        |  |  |
|                         | SEC0090   | 25ml       |  |  |
|                         | SEC0091   | 100ml      |  |  |
|                         | SEC0092   | 500ml      |  |  |
| Smartarose 4FF          | SEC0093   | 1L         |  | Smartarose 4FF and 6FF are based on highly cross-linked 4B and 6B agarose beads, exhibiting excellent physical and chromatographic properties. Due to their adaptability to various production processes and great batch repeatability, Smartarose 4FF and 6FF are widely utilized in industrial production. They are particularly effective for purifying large molecules, such as plasmids, recombinant proteins, and viruses. |
|                         | SEC0094   | 10L        |  |  |
|                         | SEC0100   | 25ml       |  |  |
|                         | SEC0101   | 100ml      |  |  |
| HiSelect Smartarose 4FF | SEC0102   | 500ml      |  |  |
|                         | SEC0103   | 1L         |  |  |
|                         | SEC0104   | 10L        |  |  |
|                         | SEC010C47 | 4.7ml,1pcs |  |  |
| HiPur Smartarose 4FF    | SEC010C20 | 20ml,1pcs  |  |  |
|                         | SEC0110   | 25ml       |  |  |
|                         | SEC0111   | 100ml      |  |  |
|                         | SEC0112   | 500ml      |  |  |
| Smartarose 6FF          | SEC0113   | 1L         |  |  |
|                         | SEC0114   | 10L        |  |  |
|                         | SEC011C47 | 4.7ml,1pcs |  |  |
|                         | SEC011C20 | 20ml,1pcs  |  |  |

### Smartdex

Smartdex is a size-exclusion chromatography medium based on dextran. Smartdex beads feature a three-dimensional network structure formed by cross-linked dextran and a cross-linking agent through ether bonds, resulting in a resin with excellent selectivity and high resolution. Beads with a higher degree of cross-linking have smaller pore sizes, allowing smaller molecules to enter the bead interiors, while larger molecules pass directly through the bead exterior.

| Products       | Separation Range(Da) (globulin) | Separation Range(Da) (linear dextran) | Swelling Volume (ml/g) | Maximum Flow Rate | Working pH |
|----------------|---------------------------------|---------------------------------------|------------------------|-------------------|------------|
| Smartdex G-10  | 700                             | 700                                   | 2-3                    | 150               | 2-13       |
| Smartdex G-15  | 1,500                           | 1,500                                 | 2.5-3.5                | 150               | 2-13       |
| Smartdex G-25  | 1,000-5,000                     | 100-5,000                             | 4-6                    | 100               | 2-13       |
| Smartdex G-50  | 1,000-30,000                    | 500-10,000                            | 9-11                   | 100               | 2-13       |
| Smartdex G-75  | 3,000-70,000                    | 1,000-50,000                          | 12-15                  | 50                | 2-10       |
| Smartdex G-100 | 4,000-150,000                   | 1,000-100,000                         | 15-20                  | 50                | 2-10       |

### Product Information

| Products               | Cat.No.  | Size | Applications  |
|------------------------|----------|------|---|
| Smartdex G-10          | SEC014S2 | 25g  | Separation and purification of polypeptide and small molecule samples.  |
|                        | SEC0141  | 100g |   |
|                        | SEC0145  | 500g |   |
|                        | SEC0146  | 1Kg  |   |
|                        | SEC0147  | 10Kg |   |
| Smartdex G-15          | SEC013S2 | 25g  |   |
|                        | SEC0131  | 100g |   |
|                        | SEC0135  | 500g |   |
|                        | SEC0136  | 1Kg  |   |
|                        | SEC0137  | 10Kg |   |
| Smartdex G-25 (Coarse) | SEC001S2 | 25g  | Smartdex G-25 is a gel filtration medium based on cross-linked dextran. Coarse particles provide a higher flow rate, while fine particles offer greater resolution. They are mainly used for separating substances of different molecular weights, such as in desalting and buffer exchange applications. |
|                        | SEC0011  | 100g |   |
|                        | SEC0015  | 500g |   |
|                        | SEC0016  | 1Kg  |   |
|                        | SEC0017  | 10Kg |   |
| Smartdex G-25 (Medium) | SEC003S2 | 25g  |   |
|                        | SEC0031  | 100g |   |
|                        | SEC0035  | 500g |   |
|                        | SEC0036  | 1Kg  |   |
|                        | SEC0037  | 10Kg |   |
| Smartdex G-25 (Fine)   | SEC005S2 | 25g  | They are mainly used for separating substances of different molecular weights, such as polypeptides, enzymes, and proteins.   |
|                        | SEC0051  | 100g |   |
|                        | SEC0055  | 500g |   |
|                        | SEC0056  | 1Kg  |   |
|                        | SEC0057  | 10Kg |   |
| Smartdex G-50 (Medium) | SEC017S2 | 25g  |   |
|                        | SEC0171  | 100g |   |
|                        | SEC0175  | 500g |   |
|                        | SEC0176  | 1Kg  |   |
|                        | SEC0177  | 10Kg |   |
| Smartdex G-50 (Fine)   | SEC021S2 | 25g  |   |
|                        | SEC02111 | 100g |   |
|                        | SEC02112 | 500g |   |
|                        | SEC02113 | 1Kg  |   |
|                        | SEC02114 | 10Kg |   |
| Smartdex G-75          | SEC012S2 | 25g  | Smartdex G-100 is commonly used for purification and separation of proteins and polysaccharides, as well as for molecular weight determination.   |
|                        | SEC0121  | 100g |   |
|                        | SEC0125  | 500g |   |
|                        | SEC0126  | 1Kg  |   |
|                        | SEC0127  | 10Kg |   |
| Smartdex G-100         | SEC015S2 | 25g  |   |
|                        | SEC0151  | 100g |   |
|                        | SEC0155  | 500g |   |
|                        | SEC0156  | 1Kg  |   |
|                        | SEC0157  | 10Kg |   |
| Smartdex LH-20         | SEC0161  | 25g  | Smartdex LH-20 can be used for separation of organic compounds with different molecular weights.  |
|                        | SEC0162  | 100g |   |
|                        | SEC0163  | 1Kg  |   |
|                        | SEC0164  | 10Kg |   |

### Agarose-dextran composite SEC resin

The Ultradex 200 Pre-packed Column is designed for high-resolution purification of large proteins. It is packed with Ultradex 200 size exclusion chromatography media, an agarose-dextran composite microsphere. This medium offers high resolution, excellent chemical stability, low non-specific binding, and easy scalability. It is ideal for efficient, high-resolution separation of monoclonal antibodies, viral vectors, and other biomolecules with molecular weights ranging from 10–600 kDa (globular proteins).

#### Example: PrePacked Column Ultradex 200 16/600




Fig. 1-19 Mouse Serum Sample Purification

### Product Information

| Products                             | Cat.No.      | Size   | Applications  |
|--------------------------------------|--------------|--------|---|
| PrePacked Column Ultradex 200 16/600 | SEC052C16/60 | 120ml  | It is ideal for efficient, high-resolution separation of monoclonal antibodies, viral vectors, and other biomolecules with molecular weights ranging from 10–600 kDa (globular proteins). |
| PrePacked Column Ultradex 200 26/600 | SEC052C26/60 | 320ml  |   |
| Ultradex 200                         | SEC0521      | 25 ml  |   |
|                                      | SEC0522      | 100 ml |   |
|                                      | SEC0523      | 500 ml |   |
|                                      | SEC0524      | 1 L    |   |
|                                      | SEC0525      | 10 L   |   |

### Desalting Column

Smartdex G-25 is effective for separating substances with molecular weights below 5000 Da, making it ideal for the rapid removal of small molecular impurities, residues, by-products, and for buffer exchange. Various desalting column formats are available to suit experimental needs, including Prepacked Desalting columns, SpinDesalt columns, and PreCap Desalting columns.




Fig. 1-20 Workflow of the Centrifugal Desalting Column.

| Products           | 0.5ml SpinDesalt Column | 3ml SpinDesalt Column | 10ml SpinDesalt Column         | P96 SpinDesalt Column                          | 8.3ml Prepacked Desalting Column | 1ml PreCap Desalting           | 5ml PreCap Desalting | 20ml HiPur Desalting            |
|--------------------|-------------------------|-----------------------|--------------------------------|--|----------------------------------|--------------------------------|----------------------|---------------------------------|
| Required Equipment | Centrifuge              |                       |                                |  | Gravity column rack              | Protein purifier               |                      |                                 |
| Loading volume     | 100–200 µl              | 350 µl–750 µl         | 1–2.5 ml                       | 60–180 µl                                      | 1.0–2.5 ml                       | 100–300 µl                     | 500–1500 µl          | 2–6 ml                          |
| Time               | 5–10 min                | 10–15 min             |                                |  | 10–20 min                        | 10–30 min                      |                      |                                 |
| Volume change      | Unchanged               |                       |                                |  | Dilute 1.5–3 times               | Dilute 1.5–2 times             |                      |                                 |
| Application        | Fast medesalting        |                       | Fast desalting of small volume | High throughput fast desalting of small volume | Fast desalting of middle volume  | Fast desalting of small volume |                      | Fast desalting of middle volume |

Fig. 1-21 Comparison of different protein desalting methods.

### Product Information

| Products                         | Cat.No.   | Size         | Applications   |
|----------------------------------|-----------|--------------|--|
| PreCap Desalting (Smartdex G-25) | SEC018C11 | 1×1ml        | They are pre-packed columns and mainly used for separating substances of different molecular weights, such as in desalting and buffer exchange applications. |
|                                  | SEC018C51 | 5×1ml        |  |
|                                  | SEC018C15 | 1×5ml        |  |
|                                  | SEC018C55 | 5×5ml        |  |
| HiSelect Desalting               | SEC003C47 | 4.7ml, 1pcs  | They are pre-packed columns that use gravity for desalting.  |
| HiPur Desalting                  | SEC003C20 | 20ml, 1pcs   |  |
| Prepacked Desalting Column       | SEC003C1  | 8.3ml, 1pcs  | They are pre-packed columns that use centrifuge for desalting.   |
|                                  | SEC003C2  | 8.3ml, 5pcs  |  |
|                                  | SEC003C3  | 8.3ml, 20pcs |  |
| SpinDesalt Column                | SEC02303  | 3T           | They are pre-packed columns that use centrifuge for desalting.   |
|                                  | SEC02301  | 50T          |  |
| SpinDesalt Column-3              | SEC04601  | 1pcs         |  |
|                                  | SEC04602  | 5pcs         |  |
|                                  | SEC04603  | 20pcs        |  |
| SpinDesalt Column-10             | SEC03201  | 1pcs         |  |
|                                  | SEC03202  | 5pcs         |  |
|                                  | SEC03203  | 20pcs        |  |
| SpinDesalt Column-P96            | SEC04701  | 1Plate       |  |
|                                  | SEC04702  | 5Plates      |  |

## 1.5 Multimodal Chromatography purification resin

Multimodal chromatography resin have more than one interaction exists between the stationary phase and the analyte. The ligands employed in these resins have different modes of interaction during purification, including size exclusion, affinity, ion exchange, and hydrophobic interactions. These diverse interactions enhance selectivity throughout the purification process.

### SmartCore 700 (Polymer matrix)

SmartCore 700 can be used in the intermediate purification and polishing of viruses and other large biomolecules. The resin features a functional core surrounded by a hydrophilic outer layer, which prevents larger molecules from entering the core. In contrast, smaller molecules and impurities can access the core, where they bind to the hydrophobic and positively charged octylamine ligands. The molecular size cut-off for proteins is approximately 690 kDa.

#### Product Information

| Products               | Cat.No.   | Size       | Applications   |
|------------------------|-----------|------------|--|
| SmartCore 700          | SEC0220   | 25ml       | SmartCore 700 multimodal chromatography resins is designed for intermediate purification and polishing of viruses and other large biomolecules. Known for its higher flow rate, SmartCore 700 could significantly improve production efficiency, and is suitable for large-scale production. |
|                        | SEC0221   | 100ml      |  |
|                        | SEC0222   | 500ml      |  |
|                        | SEC0223   | 1L         |  |
|                        | SEC0224   | 10L        |  |
| PreCap SmartCore 700   | SEC022C11 | 1×1ml      |  |
|                        | SEC022C51 | 5×1ml      |  |
|                        | SEC022C15 | 1×5ml      |  |
|                        | SEC022C55 | 5×5ml      |  |
| HiSelect SmartCore 700 | SEC022C47 | 4.7ml,1pcs |  |
| HiPur SmartCore 700    | SEC022C20 | 20ml,1pcs  |  |

### Smac Core 700 (Highly rigid agarose spheres matrix)

Smac Core 700 can be used in the intermediate purification and polishing of viruses and other large biomolecules. The resin features a functional core surrounded by a hydrophilic outer layer, which prevents larger molecules from entering the core. In contrast, smaller molecules and impurities can access the core, where they bind to the hydrophobic and positively charged octylamine ligands. The molecular size cut-off for proteins is approximately 690 kDa.

#### Product Information

| Products             | Cat.No.   | Size  | Applications   |
|----------------------|-----------|-------|--|
| Smac Core 700        | SEC0281   | 25ml  | Smac Core 700 utilizes a core bead design, with the core activated by octylamine and a ligand-free outer layer. This multimodal chromatography medium is ideal for intermediate purification and polishing of large biomolecules in flow-through mode. |
|                      | SEC0282   | 100ml |  |
|                      | SEC0283   | 500ml |  |
|                      | SEC0284   | 1L    |  |
| PreCap Smac Core 700 | SEC028C11 | 1×1ml |  |
|                      | SEC028C51 | 5×1ml |  |
|                      | SEC028C15 | 1×5ml |  |
|                      | SEC028C55 | 5×5ml |  |

### Smac Core 400 (Highly rigid agarose spheres matrix)

Smac Core 400 can be used in the intermediate purification and polishing of viruses and other large biomolecules. The resin features a functional core surrounded by a hydrophilic outer layer, which prevents larger molecules from entering the core. In contrast, smaller molecules and impurities can access the core, where they bind to the hydrophobic and positively charged octylamine ligands. The molecular size cut-off for proteins is approximately 400 kDa.

#### Product Information

| Products               | Cat.No.   | Size      | Applications  |
|------------------------|-----------|-----------|---|
| Smac Core 400          | SEC0291   | 25ml      | Designed for the purification of albumin, interferons, nucleotide-requiring enzymes, α2-macroglobulin, coagulation factors, and other related biomolecules. |
|                        | SEC0292   | 100ml     |   |
|                        | SEC0293   | 500ml     |   |
|                        | SEC0294   | 1L        |   |
|                        | SEC0295   | 10L       |   |
| PreCap Smac Core 400   | SEC029C11 | 1×1ml     |   |
|                        | SEC029C51 | 5×1ml     |   |
|                        | SEC029C15 | 1×5ml     |   |
|                        | SEC029C55 | 5×5ml     |   |
| HiSelect Smac Core 400 | SEC029C47 | 4.7ml     |   |
| HiPur Smac Core 400    | SEC029C20 | 20ml,1pcs |   |

### Smac MMA

Smac MMA is a multimodal strong anion exchange resin with ionic interaction, hydrogen bond, hydrophobic interaction. Smac MMA can remove key impurities such as DNA, host cell proteins (HCP), leached protein A, aggregates, and viruses in post-protein A purification steps. It can also be used for fine purification of other biomolecules, such as the removal of dimers, polymers, host cell proteins, and nucleic acids.

#### Product Information

| Products          | Cat.No.  | Size       | Applications  |
|-------------------|----------|------------|---|
| Smac MMA          | SI031025 | 25ml       | Smac MMA can be used for moderate and fine purification of monoclonal antibodies, as well as for the fine purification of other biomolecules. |
|                   | SI031100 | 100ml      |   |
|                   | SI031500 | 500ml      |   |
|                   | SI03101L | 1L         |   |
|                   | SI03110L | 10L        |   |
| lexCap Smac MMA   | SI031C11 | 1×1ml      |   |
|                   | SI031C51 | 5×1ml      |   |
|                   | SI031C15 | 1×5ml      |   |
|                   | SI031C55 | 5×5ml      |   |
| HiSelect Smac MMA | SI031C47 | 4.7ml,1pcs |   |
| HiPur Smac MMA    | SI031C20 | 20ml,1pcs  |   |

## Smac MMC

Smac MMC is a multimodal weak cation exchange resin with ionic interaction, hydrogen bond, hydrophobic interaction. It can be used for fast and efficient protein purification. The Smac MMC has the following characteristics:

Unique selectivity compared to conventional ion exchange media.

High dynamic binding capacity under high conductivity conditions.



Easy to optimize and scale up.

Designed to meet industrial demands for supply security, consistent performance, and regulatory support.

### Product Information

| Products             | Cat.No.  | Size       | Applications   |
|----------------------|----------|------------|--|
| Smac MMC             | SI033025 | 25ml       | Smac MMC can be used for fast and efficient protein purification.  |
|                      | SI033100 | 100ml      |  |
|                      | SI033500 | 500ml      |  |
|                      | SI03301L | 1L         |  |
|                      | SI03310L | 10L        |  |
| lexCap Smac MMC      | SI033C11 | 1×1ml      |  |
|                      | SI033C51 | 5×1ml      |  |
|                      | SI033C15 | 1×5ml      |  |
|                      | SI033C55 | 5×5ml      |  |
| HiSelect Smac MMC    | SI033C47 | 4.7ml,1pcs |  |
| HiPur Smac MMC       | SI033C20 | 20ml,1pcs  |  |
| Smac MMC 40          | SI032025 | 25ml       | Smac MMC 40 can be used for moderate and fine purification of monoclonal antibodies. The product has high binding ability selectivity over a wide range of pH and salt concentrations. |
|                      | SI032100 | 100ml      |  |
|                      | SI032500 | 500ml      |  |
|                      | SI03201L | 1L         |  |
|                      | SI03210L | 10L        |  |
| lexcap Smac MMC 40   | SI032C11 | 1×1ml      |  |
|                      | SI032C15 | 1×5ml      |  |
|                      | SI032C51 | 5×1ml      |  |
|                      | SI032C55 | 5×5ml      |  |
| HiSelect Smac MMC 40 | SI032C47 | 4.7ml,1pcs |  |
| HiPur Smac MMC 40    | SI032C20 | 20ml,1pcs  |  |

# Polymer-Based Chromatography Resin

Phmac series hydrophobic interaction chromatography (HIC)resin

Phmac series ion exchange chromatography (IEC) resin

Smart PS series ion exchange chromatography (IEC)resin

Smart PS series reversed phase chromatography (RPC)resin

# 02

Polymer chromatography resins are primarily used in the preparation and purification of small molecules such as peptides, chemical drugs, and natural products, as well as biomacromolecules. The success of production depends on economical and efficient purification and separation solutions. Resin and process selection must ensure safety and effectiveness, enabling the efficient removal of process-related impurities, including host cell proteins, residual DNA, viruses, and endotoxins, while preserving the activity and recovery yield of the target product. The polymer chromatography resins developed by smart-lifesciences incorporate technologies in matrix synthesis, surface hydrophilicity modification, and ligand coupling, conferring excellent hydrophilicity, high binding capacity, long service life, and high flow performance. These resins significantly improve downstream purification efficiency, reduce overall process costs, and deliver substantial economic value.

## 2.1 Phmac Series HIC Resin

Smart-lifesciences has launched a product featuring monodisperse polymethacrylate microspheres as the base matrix, coupled with hydrophobic ligands such as butyl and phenyl groups. This product operates by adsorbing biomolecules under high-salt conditions and eluting them under low-salt conditions. It exhibits excellent biocompatibility, high binding capacity, strong chemical stability, and superior mechanical strength. It can be widely applied in the separation and purification of antibodies, recombinant proteins, vaccines, peptides, nucleic acids, and more.

### Product Advantages

1. Monodisperse microspheres offer superior batch-to-batch consistency.
2. High mechanical strength supports higher flow rates and operating pressures, improving purification efficiency.
3. Low compression coefficient ensures a stable column bed.
4. Excellent alkali resistance extends resin lifespan and reduces production costs.
5. Comprehensive product specifications, supported by professional customization services, meet diverse application requirements.

### Product Description

| Product   | Butyl Phmac 30 Beads  | Phenyl Phmac 30 Beads    | Butyl Phmac 60 Beads             | Phenyl Phmac 60 Beads            |
|---|---|--------------------------|----------------------------------|----------------------------------|
| Material  | PGMA  |                          |                                  |                                  |
| Type  | Butyl   | Phenyl                   | Butyl                            | Phenyl                           |
| Particle Size (µm)                                      | ~30µm   |                          | ~60µm                            |                                  |
| Pore Size (nm) <sup>1</sup>                             | 50(XS)/100(S)/200(M)/500(L)nm   |                          |                                  |                                  |
| Recommended Flow Rate Range                             | 50-300cm/h  |                          | 100-600 cm/h                     |                                  |
| Static Binding Capacity <sup>2</sup> (mg lys/ml medium) | ~ 45(XS) ~ 30(S) ~ 20(M)  | ~ 55(XS) ~ 30(S) ~ 20(M) | ~ 40(XS) ~ 25(S) ~ 15(M) ~ 10(L) | ~ 50(XS) ~ 25(S) ~ 15(M) ~ 10(L) |
| pH Stability <sup>2</sup>                               | 2-13  |                          |                                  |                                  |
| Maxi Pressure   | ≤1 Mpa  |                          | ≤0.8 Mpa                         |                                  |
| Chemical Stability                                      | Compatible with water and mixed solvents containing ethanol, acetonitrile, acetone, etc. Commonly used buffer systems: Tris, phosphate, acetate buffer systems. |                          |                                  |                                  |
| Application Range                                       | Widely applied in the separation and purification of antibodies, recombinant proteins, vaccines, peptides, nucleic acids, etc.                                  |                          |                                  |                                  |
| Storage Conditions                                      | 4-35 °C, in°C 50 % (v/v) 20 % ethanol.  |                          |                                  |                                  |

1. Dynamic Binding Capacity is measured using Lysozyme, with data for different pore sizes provided in parentheses.  
2. pH Stability Range refers to the applicable pH intervals for operation, regeneration, and Cleaning In Place (CIP).

## Application Case 1: Capture of Fusion Protein Supernatant Expressed in E. coli – Butyl Phmac Beads




Fig. 2-1 Capture of Fusion Protein Supernatant Expressed in E. coli – Butyl Phmac Beads

### Conclusion:

- Post-capture sample purity is significantly improved.
- Two analytical methods indicate/demonstrate that purity is increased to greater than 92%.
- The yield exceeds 90%, meeting the expected targets/specifications for both yield and purity.

## Application Case 2: Purification of Fusion Protein after Cleavage – Butyl Phmac Beads




Fig. 2-2 Purification of Fusion Protein after Cleavage – Butyl Phmac Beads

### Conclusion:

- No target protein was detected in the flow-through fraction during sample loading.
- Purity > 90%, yield > 90%.
- CIP (Cleaning In Place) indicates almost no loss of the main component.

Product Information

| Products                  | Product Code | Resin Specifications        | Pre-packed columns | Application  |
|---------------------------|--------------|-----------------------------|--------------------|--|
| Butyl Phmac 30 Beads(XS)  | SH054        | 5ml,25ml,100ml,500ml,1L,10L | 1ml,5ml,4.7ml,20ml | Widely used in the separation and purification of antibodies, recombinant proteins, plasmids, viruses, vaccines, peptides, nucleic acids, etc. |
| Butyl Phmac 30 Beads(S)   | SH051        |                             |                    |  |
| Phenyl Phmac 30 Beads(XS) | SH055        |                             |                    |  |
| Phenyl Phmac 30 Beads(S)  | SH052        |                             |                    |  |
| Butyl Phmac Beads         | SH019        |                             |                    |  |
| Butyl Phmac 60 Beads(S)   | SH040        |                             |                    |  |
| Phenyl Phmac 60 Beads(XS) | SH039        |                             |                    |  |
| Phenyl Phmac 60 Beads(S)  | SH049        |                             |                    |  |

## 2.2 Phmac Series Ion Exchange Chromatography (IEC) Resin

Phmac Ion Exchange Resins employ a hydrophilic polymethacrylate matrix with a hydroxyl-rich surface, enabling efficient coupling with various ion exchange functional groups, including DEAE, Q, CM, and SP. These resins combine high mechanical strength, excellent hydrophilicity, and extremely low nonspecific adsorption, supporting rapid mass transfer and high binding capacity even under high flow-rate conditions. Designed specifically for the fine purification of biomacromolecules (e.g., antibodies, vaccines, recombinant proteins) and active components from traditional Chinese medicine, they maintain outstanding stability under stringent cleaning and regeneration conditions, offering an ideal chromatography solution for enhancing downstream purification efficiency and product purity.

**Key Features of the Phmac Series Ion Exchange Resins:**

- High dynamic binding capacity. Two particle size options (30 μm and 60 μm) to meet different chromatography requirements.
- Low compression coefficient ensures a stable column bed.
- High rigidity, enabling operation at higher flow rates and pressures to enhance process efficiency.
- Excellent alkali resistance extends resin lifespan and reduces production costs.

Product Description

| Resin                      | Phmac 30 SP (S)/(M)/(L)   | Phmac 30 CM (S)/(M)/(L) | Phmac 30 Q (S)/(M)/(L)                          | Phmac 30 DEAE (S)/(M)/(L)                        | Phmac 60 SP (XS)/(S)/(M)/(L)       | Phmac 60 CM (XS)/(S)/(M)/(L) | Phmac 60 Q (XS)/(S)/(M)/(L)                     | Phmac 60 DEAE (XS)/(S)/(M)/(L)                   |
|----------------------------|---|-------------------------|---|--|------------------------------------|------------------------------|---|--|
| Ion Exchange Type          | Strong Cation Exchange  | Weak Cation Exchange    | Strong Anion Exchange                           | Weak Anion Exchange                              | Strong Cation Exchange             | Weak Cation Exchange         | Strong Anion Exchange                           | Weak Anion Exchange                              |
| Ligand                     | -SO <sub>3</sub> <sup>-</sup>   | -COO <sup>-</sup>       | -N <sup>+</sup> (CH <sub>3</sub> ) <sub>3</sub> | -N <sup>+</sup> H(CH <sub>3</sub> ) <sub>2</sub> | -SO <sub>3</sub> <sup>-</sup>      | -COO <sup>-</sup>            | -N <sup>+</sup> (CH <sub>3</sub> ) <sub>3</sub> | -N <sup>+</sup> H(CH <sub>3</sub> ) <sub>2</sub> |
| Matrix                     | PGMA  |                         |   |  | PGMA                               |                              |   |  |
| Particle Size              | ~30μm   |                         |   |  | ~60μm                              |                              |   |  |
| Pore Size                  | 100(S)/200(M)/500(L) nm   |                         |   |  | 50(XS)/100(S)/200(M)/500(L) nm     |                              |   |  |
| Ion Capacity               | 0.15-0.25 mmol/ml   |                         |   |  |                                    |                              |   |  |
| DBC                        | S,M≥80 mg/ml,L≥50 mg/ml (lysozyme)  |                         | S,M≥80 mg/ml,L≥50 mg/ml (BSA)                   |  | S,M≥80 mg/ml,L≥50 mg/ml (lysozyme) |                              | S,M≥80 mg/ml,L≥50 mg/ml (BSA)                   |  |
| Maximum Pressure Tolerance | 1MPa  |                         |   |  | 0.5MPa                             |                              |   |  |
| CIP                        | 0.5M NaOH   |                         |   |  |                                    |                              |   |  |
| Recommended Flow Rate      | 50-300cm/h  |                         |   |  | 150-750cm/h                        |                              |   |  |
| pH Stability               | 2-12  |                         |   |  |                                    |                              |   |  |
| Chemical Stability         | Compatible with all common buffers, 1M acetic acid, 1M sodium hydroxide, 1M hydrochloric acid, 70% ethanol, 30% isopropanol, 30% acetonitrile, 1% SDS, 6M guanidine hydrochloride, 8M urea, and other common organic solvents. Avoid exposure to strong oxidizing agents. |                         |   |  |                                    |                              |   |  |
| Storage Buffer             | 20% ethanol   |                         |   |  |                                    |                              |   |  |
| Storage Temperature        | 4-30°C  |                         |   |  |                                    |                              |   |  |




Fig. 2-3 Purification of Captured Plasmid Samples - Butyl Phmac 60 Beads (M)




Fig. 2-4 Agarose Gel Electrophoresis of Purified Plasmid Samples

**Conclusion:** Purity of the sample was significantly improved after polishing, with scDNA purity reaching above 95%.

Product Information

| Products        | Product Code | Resin Specifications         | Pre-packed columns | Applications  |
|-----------------|--------------|------------------------------|--------------------|---|
| Phmac 30 SP(S)  | SI087        | 5ml,25ml,100ml,500ml, 1L,10L | 1ml,5ml,4.7ml,20ml | Polymethacrylate-Based Cation Exchange Resin: Available in 30μm particle size with three pore size options: 100 nm (S), 200 nm (M), and 500 nm (L); and in 60μm particle size with four pore size options: 50 nm (XS), 100 nm (S), 200 nm (M), and 500 nm (L). Suitable for industrial-scale separation and purification of biomacromolecules such as antibodies, recombinant proteins, and blood products, and also applicable for removing charge variants and host cell proteins from protein samples. |
| Phmac 30 SP(M)  | SI092        |                              |                    |   |
| Phmac 30 SP(L)  | SI096        |                              |                    |   |
| Phmac 60 SP(XS) | SI123        |                              |                    |   |
| Phmac 60 SP(S)  | SI085        |                              |                    |   |
| Phmac 60 SP(M)  | SI105        |                              |                    |   |
| Phmac 60 SP(L)  | SI109        |                              |                    |   |
| Phmac 30 CM(S)  | SI089        |                              |                    |   |
| Phmac 30 CM(M)  | SI093        |                              |                    |   |
| Phmac 30 CM(L)  | SI097        |                              |                    |   |
| Phmac 60 CM(XS) | SI122        |                              |                    |   |
| Phmac 60 CM(S)  | SI086        |                              |                    |   |
| Phmac 60 CM(M)  | SI106        |                              |                    |   |
| Phmac 60 CM(L)  | SI110        |                              |                    |   |

| Products          | Product Code | Resin Specifications        | Pre-packed columns | Applications   |
|-------------------|--------------|-----------------------------|--------------------|--|
| Phmac 30 Q(S)     | SI081        | 5ml,25ml,100ml,500ml,1L,10L | 1ml,5ml,4.7ml,20ml | Polymethacrylate-Based Anion Exchange Resin: Available in 30µm particle size with three pore size options: 100 nm (S), 200 nm (M), and 500 nm (L); and in 60µm particle size with four pore size options: 50 nm (XS), 100 nm (S), 200 nm (M), and 500 nm (L). Suitable for industrial-scale separation and purification of biomacromolecules such as antibodies, vaccines, virus-like particles, recombinant proteins, and blood products. Also applicable for nucleic acid purification and removal of host nucleic acids and host proteins from protein samples. |
| Phmac 30 Q(M)     | SI091        |                             |                    |  |
| Phmac 30 Q(L)     | SI095        |                             |                    |  |
| Phmac 60 Q(XS)    | SI124        |                             |                    |  |
| Phmac 60 Q(S)     | SI071        |                             |                    |  |
| Phmac 60 Q(M)     | SI104        |                             |                    |  |
| Phmac 60 Q(L)     | SI108        |                             |                    |  |
| Phmac 30 DEAE(S)  | SI080        |                             |                    |  |
| Phmac 30 DEAE(M)  | SI090        |                             |                    |  |
| Phmac 30 DEAE(L)  | SI094        |                             |                    |  |
| Phmac 60 DEAE(XS) | SI102        |                             |                    |  |
| Phmac 60 DEAE(S)  | SI072        |                             |                    |  |
| Phmac 60 DEAE(M)  | SI103        |                             |                    |  |
| Phmac 60 DEAE(L)  | SI107        |                             |                    |  |

### 2.3

## Smart PS Series Ion Exchange Chromatography (IEC) Resin

Smart PS Series Ion Exchange Resins are based on a polystyrene-divinylbenzene copolymer matrix, which undergoes excellent hydrophilic modification before coupling with ion-exchange functional groups. These resins combine high mechanical strength with extremely low nonspecific adsorption, enabling rapid mass transfer and high binding capacity even at high flow rates. They significantly enhance purification efficiency and are suitable for high-resolution purification of biomolecules such as plasmids, peptides, nucleotides, and insulin. With outstanding stability under harsh cleaning and regeneration conditions, they serve as an ideal chromatography solution to improve downstream purification efficiency and product purity.

#### Advantages of the Smart PS Series Ion Exchange Resins:

- Monodisperse polystyrene-divinylbenzene matrix provides high mechanical strength and high resolution.
- Hydrophilic surface modification ensures low nonspecific adsorption.
- Low compression coefficient maintains a stable column bed.
- Low backpressure supports high-flow-rate rapid purification.
- Excellent chemical resistance extends resin service life.

### Product Description

| Product  | Smart PS 15Q   | Smart PS 15S       | Smart PS 30Q                                    | Smart PS 30S       |
|--|--|--------------------|---|--------------------|
| Material   | PS/DVB   |                    |   |                    |
| Type   | -N <sup>+</sup> (CH <sub>3</sub> ) <sub>3</sub>  | -SO <sub>3</sub> H | -N <sup>+</sup> (CH <sub>3</sub> ) <sub>3</sub> | -SO <sub>3</sub> H |
| Particle Size  | ~15µm  |                    | ~30µm   |                    |
| Recommended Flow Rate Range                          | 150-750cm/h  |                    | 200-950cm/h                                     |                    |
| Dynamic Binding Capacity <sup>1</sup> (mg/ml medium) | ~55mg/ml(BSA)  | ~70mg/ml(Lysozyme) | ~45mg/ml(BSA)                                   | ~60mg/ml(Lysozyme) |
| pH Stability <sup>2</sup>                            | 2-13   |                    |   |                    |
| Pressure Tolerance                                   | ≤6Mpa  |                    | ≤3Mpa   |                    |
| Chemical Stability                                   | Compatibility: Compatible with water and mixed solvents containing ethanol, acetonitrile, acetone, etc. Commonly Used Buffer Systems: Tris, phosphate, acetate buffer systems. |                    |   |                    |
| Application Range                                    | Widely applied in the separation and purification of antibodies, recombinant proteins, vaccines, peptides, nucleic acids, insulin, etc.  |                    |   |                    |
| Storage Conditions                                   | 4 °C to 35 °C, in 50 % (v/v) 20 % ethanol.   |                    |   |                    |

1. Dynamic Binding Capacity is measured using Lysozyme, with data for different pore sizes provided in parentheses.

2. pH Stability Range refers to the applicable pH intervals for operation, regeneration, and Cleaning In Place (CIP).

## Application Case : Ion Exchange Purification of GLP-1 after Side Chain Modification




Fig. 2-5 Ion Exchange Purification of GLP-1 after Side Chain Modification

Purification using an alkaline buffer system can effectively improve the purity of the crude product after structural modification:

| Products     | Loading Capacity (based on pure peptide weight) | Crude Product Purity | Collected Main Peak Purity | Yield |
|--------------|---|----------------------|----------------------------|-------|
| Smart PS 15Q | 20mg/ml   | 59%                  | 93.8%                      | 96%   |

### Product Information

| Products     | Product Code | Resin Specifications        | Pre-packed columns | Applications   |
|--------------|--------------|-----------------------------|--------------------|--|
| Smart PS 15Q | SI066        | 5ml,25ml,100ml,500ml,1L,10L | 1ml,5ml,4.7ml,20ml | Widely used in the separation and purification of antibodies, recombinant proteins, vaccines, peptides, nucleic acids, insulin, etc. |
| Smart PS 15S | SI068        |                             |                    |  |
| Smart PS 30Q | SI069        |                             |                    |  |
| Smart PS 30S | SI070        |                             |                    |  |

Medium-/High-Pressure Stainless-Steel Analytical/Preparative Columns: 4.6/10/20 × 150/250 mm

### 2.4

## Smart PS Series RPC Resin

Smart PS Series Reversed-Phase Polymeric Resins are monodisperse resins based on a highly cross-linked PS/DVB (polystyrene-divinylbenzene) matrix. They are available in multiple particle sizes and pore sizes, making them suitable for various chromatography purification fields and different purification stages. Due to their abundant phenyl structures, the Smart PS series exhibits strong hydrophobicity and excellent reversed-phase retention capability. The highly cross-linked PS/DVB structure confers outstanding mechanical strength and chemical stability, resulting in superior pressure resistance and compatibility with medium-to high-pressure DAC (Dynamic Axial Compression) column packing. In addition, these resins tolerate extreme acidic and alkaline conditions as well as regeneration with mixed organic solvents. They are particularly suitable for the separation and purification of highly contaminated samples, such as fermentation broths, and for samples requiring repeated regeneration under extreme conditions.

#### Advantages of the Smart PS Series Reversed-Phase Resins:

- Monodisperse PS/DVB matrix ensures high mechanical strength and high resolution.
- Low compression coefficient maintains a stable column bed.
- Low backpressure supports high-flow-rate rapid purification.
- Excellent chemical resistance extends resin service life.

Product Information

| Product                     | Smart PS 10-300   | Smart PS 10-500 | Smart PS 15-300 | Smart PS 15-500 |
|-----------------------------|---|-----------------|-----------------|-----------------|
| Material                    | PS/DVB  |                 |                 |                 |
| Type                        | Phenyl  |                 |                 |                 |
| Particle Size               | ~10µm   |                 | ~15µm           |                 |
| Pore Size (Å)               | 300   | 500             | 300             | 500             |
| Recommended Flow Rate Range | 150-750cm/h   |                 | 200-950cm/h     |                 |
| pH Stability <sup>2</sup>   | 2-13  |                 |                 |                 |
| Swelling Factor             | ≤3% Methanol  |                 | ≤5% Methanol    |                 |
| Pressure Tolerance          | ≤10Mpa  |                 |                 |                 |
| Chemical Stability          | Compatibility: Compatible with water and mixed solvents containing ethanol, acetonitrile, acetone, etc. Commonly Used Buffer Systems: Tris, phosphate, acetate buffer systems.  |                 |                 |                 |
| Regeneration                | Column Cleaning with Solvents of Increasing Non-Polarity: For example, ethanol/methanol/acetonitrile-water mixtures, followed by pure ethanol/methanol/acetonitrile. Alkaline Cleaning: it is recommended to rinse with 0.1-0.5 M NaOH/60 % methanol solution for 3-5 column volumes (CV), then rinse with 0.1-0.5 M HCl or acetic acid/60 % methanol solution for 3-5 CV, and finally equilibrate with the mobile phase. This is because protein aggregates and other contaminants are readily disrupted and eluted in acidic or alkaline solutions. |                 |                 |                 |
| Application Range           | Widely used in the separation and purification of antibodies, recombinant proteins, vaccines, peptides, nucleic acids, insulin, etc.  |                 |                 |                 |
| Storage Conditions          | 4 °C to 35 °C, in 50 % (v/v) 20 % ethanol.  |                 |                 |                 |

Note: Due to the strong hydrophobicity of polymer-based reversed-phase resins, operation with a purely aqueous mobile phase is not recommended. During chromatographic separation, ensure that the column remains continuously wetted and that the mobile phase does not run dry, as drying may cause cracking of the column bed. Such cracks can lead to uneven distribution of the sample and mobile phase upon rewetting, resulting in poor reproducibility and compromised purification performance.

Application Case: Ion Exchange Purification of GLP-1 after Side Chain Modification




Fig. 2-6 Ion Exchange Purification of GLP-1 after Side Chain Modification

Purification using an alkaline buffer/acetonitrile system can effectively enhance the purity of the crude product following structural modification.

| Fraction         | Peak Area | Dilution Factor | Loading/Collection Volume (ml) | Chromatographic Purity | Yield  |
|------------------|-----------|-----------------|--------------------------------|------------------------|--------|
| Original Sample  | 3360      | 5               | 5                              | 72%                    | NA     |
| Collected Sample | 8992.8    | 3               | 3                              | 97.69%                 | 96.35% |

Product Information

| Products        | Product Code | Resin Specifications         | Applications   |
|-----------------|--------------|------------------------------|--|
| Smart PS 10-300 | SH037        | 5ml,25ml,100ml,500ml,1L, 10L | Reverse-phase chromatographic separation of small molecule compounds, peptides, low-molecular-weight proteins, and other biomolecules. |
| Smart PS 10-500 | SH038        |                              |  |
| Smart PS 15-300 | SH035        |                              |  |
| Smart PS 15-500 | SH036        |                              |  |

Stainless-Steel Analytical/Preparative Columns: 4.6/10/20 × 150/250 mm

# Other Chromatography Resin

Smart C8/C18

MonosilicaC8/C18

Hydroxyapatite (CHT)

# 03

In response to the growing demand for spherical silica, we have developed an optimized large-scale production process that provides a high level of process control and can be applied across all smart-lifesciences silica products without compromising quality. Key parameters, including particle morphology, pore size, particle size distribution, silica purity, and surface chemistry, directly influence chromatographic performance. Therefore, to develop an excellent process, it is necessary to evaluate and optimize all these parameters to ensure batch-to-batch reproducibility of the product.

### 3.1 Smart C8/C18

Smart C8/C18 is a high-performance functionalized spherical silica specifically designed to meet the stringent demands of peptide purification. Engineered for both analytical and preparative chromatography, it offers exceptional resilience, efficiency, reproducibility, and scalability, making it an ideal solution for pharmaceutical, biotechnology, and research applications.


| Products         | Particle Size (µm) | Pore Size (Å) | Specific Surface Area (m <sup>2</sup> /g) | Carbon Content (%) | Bonded Phase | pH    |
|------------------|--------------------|---------------|---|--------------------|--------------|-------|
| Smart C8 10-100  | 10                 | 100           | 360-440                                   | 14                 | C8           | 2-9.5 |
| Smart C18 10-100 | 10                 | 100           | 360-440                                   | 22                 | C18          | 2-9.5 |
| Smart C8 10-120  | 10                 | 120           | 290-340                                   | 12                 | C8           | 2-9.5 |
| Smart C18 10-120 | 10                 | 120           | 290-340                                   | 20                 | C18          | 2-9.5 |
| Smart C8 10-150  | 10                 | 150           | 220-270                                   | 10                 | C8           | 2-9.5 |
| Smart C18 10-150 | 10                 | 150           | 220-270                                   | 18                 | C18          | 2-9.5 |

#### Application Case

##### Step 1 Purification

Resin: Smart 18 10-100, acetic acid-acetonitrile

Preparative Chromatogram




| Sample          | RRT0.89 | RRT0.92 | Pooled Fraction Purity | RRT1.06 | Yield |
|-----------------|---------|---------|------------------------|---------|-------|
| Standard        | /       | /       | 99.288                 | /       | /     |
| Crude Product   | 0.302   | 0.142   | 59.529                 | 0.234   | /     |
| Pooled Fraction | 0.098   | 0.022   | 98.861                 | 0.295   | 93.8% |

##### Step 2 Purification

Resin: Smart 18 10-100, ammonium acetate-acetonitrile gradient

Preparative Chromatogram



| Sample           | RRT0.89 | RRT0.92 | Pooled Fraction Purity | RRT1.05 | Yield |
|------------------|---------|---------|------------------------|---------|-------|
| Standard         | /       | /       | 99.288                 | /       | /     |
| Loading Solution | 0.135   | 0.032   | 98.861                 | 0.244   | /     |
| Pooled Fraction  | 0.07    | 0.019   | 99.68                  | 0.111   | 94.7% |

### Product Information

| Products         | Product Code | Resin Specifications | Pre-Packed Column Specifications |
|------------------|--------------|----------------------|----------------------------------|
| Smart C18 10-100 | SH020        | 25g,100g,500g,1kg    | 4.6/10/20mm×250mm                |
| Smart C8 10-100  | SH021        |                      |                                  |
| Smart C8 10-150  | SH030        |                      |                                  |
| Smart C18 10-150 | SH031        |                      |                                  |
| Smart C8 10-120  | SH043        |                      |                                  |
| Smart C18 10-120 | SH044        |                      |                                  |

Custom pre-packed column sizes (length and diameter adjustable) available upon request.

### 3.2 MonoSilica C8/C18

MonoSilica C8/C18 is a versatile, high-purity silica-based chromatography resin designed for the purification of small molecules, short nucleic acids, and peptides. It delivers well-balanced performance across key parameters, including separation efficiency, loading capacity, recovery, and mechanical strength.

| Products              | Bonded Phase | Particle Size (µm) | Pore Size (Å) | Product Code |
|-----------------------|--------------|--------------------|---------------|--------------|
| MonoSilica C8 10-300  | C8           | 10                 | 300           | SH022        |
| MonoSilica C18 10-300 | C18          | 10                 | 300           | SH023        |
| MonoSilica C8 10-100  | C8           | 10                 | 100           | SH024        |
| MonoSilica C18 10-100 | C18          | 10                 | 100           | SH025        |
| Smart AQC8 10-100     | C8           | 10                 | 100           | SH028        |
| Smart AQC18 10-100    | C18          | 10                 | 100           | SH029        |

### 3.3 Hydroxyapatite( CHT)

Hydroxyapatite(CHT I , CHT II) is a unique, rigid composite chromatography medium for downstream purification that combines ion-exchange and metal-affinity mechanisms, providing separation capabilities not achievable with conventional media. The medium can capture target proteins by calcium ions or phosphate-associated positively charged groups. It is manufactured from inorganic materials sintered at high temperatures, Hydroxyapatite contains no plant- or animal-derived components. It is widely applicable to the purification of a broad range of biomolecules, including monoclonal and polyclonal antibodies, bispecific antibodies and antibody fragments, recombinant proteins, vaccines, enzymes, nucleic acids, viral particles, and other biologics. Hydroxyapatite is available in two types—CHT I and CHT II, each offered in three particle size options. CHT I is well suited for the purification of acidic proteins, while CHT II, with its larger pore structure, is optimized for large biomolecules (>400 kDa) such as IgM, viruses, virus-like particles (VLPs), and plasmids. Hydroxyapatite media are versatile and can be applied across all stages of purification, from initial capture to high-resolution polishing.

Product Information

| Products                      | Cat.No.  | Size  | Applications  |
|-------------------------------|----------|-------|---|
| Hydroxyapatite 60( CHT I 60)  | SA111005 | 5g    | It can be applied to the purification of various biomacromolecules, including monoclonal and polyclonal antibodies, bispecific antibodies, antibody fragments, recombinant proteins, vaccines, enzymes, and nucleic acids. It is particularly suitable for the purification of acidic proteins. |
|                               | SA111020 | 20g   |   |
|                               | SA111100 | 100g  |   |
|                               | SA111500 | 500g  |   |
|                               | SA1111KG | 1Kg   |   |
| Hydroxyapatite 40( CHT I 40)  | SA106005 | 5g    |   |
|                               | SA106020 | 20g   |   |
|                               | SA106100 | 100g  |   |
|                               | SA106500 | 500g  |   |
|                               | SA1061KG | 1Kg   |   |
| Hydroxyapatite 20( CHT I 20)  | SA113005 | 5g    |   |
|                               | SA113020 | 20g   |   |
|                               | SA113100 | 100g  |   |
|                               | SA113500 | 500g  |   |
|                               | SA1131KG | 1Kg   |   |
| Hydroxyapatite 60( CHT II 60) | SA112005 | 5g    | CHT II resins feature larger pore sizes, making them suitable for purifying large-particle biomolecules (>400 kDa), such as IgM, viruses, VLPs, and plasmids.   |
|                               | SA112020 | 20g   |   |
|                               | SA112100 | 100g  |   |
|                               | SA112500 | 500g  |   |
|                               | SA1121KG | 1Kg   |   |
| Hydroxyapatite 40( CHT II 40) | SA107005 | 5g    |   |
|                               | SA107020 | 20g   |   |
|                               | SA107100 | 100g  |   |
|                               | SA107500 | 500g  |   |
|                               | SA1071KG | 1Kg   |   |
| Hydroxyapatite 20( CHT II 20) | SA114005 | 5g    |   |
|                               | SA114020 | 20g   |   |
|                               | SA114100 | 100g  |   |
|                               | SA114500 | 500g  |   |
|                               | SA1141KG | 1Kg   |   |
| PreCap CHT II 20              | SA114C11 | 1*1ml | The prepacked columns are supplied with standard fittings and are compatible with a wide range of medium-pressure chromatography systems, including ÄKTA™, ensuring easy and reliable operation.  |
|                               | SA114C15 | 1*5ml |   |
|                               | SA114C51 | 5*1ml |   |
|                               | SA114C55 | 5*5ml |   |
| PreCap CHT II 40              | SA107C11 | 1*1ml |   |
|                               | SA107C15 | 1*5ml |   |
|                               | SA107C51 | 5*1ml |   |
|                               | SA107C55 | 5*5ml |   |
| PreCap CHT II 60              | SA112C11 | 1*1ml |   |
|                               | SA112C15 | 1*5ml |   |
|                               | SA112C51 | 5*1ml |   |
|                               | SA112C55 | 5*5ml |   |
| PreCap CHT I 20               | SA113C11 | 1*1ml |   |
|                               | SA113C15 | 1*5ml |   |
|                               | SA113C51 | 5*1ml |   |
|                               | SA113C55 | 5*5ml |   |
| PreCap CHT I 40               | SA106C11 | 1*1ml |   |
|                               | SA106C15 | 1*5ml |   |
|                               | SA106C51 | 5*1ml |   |
|                               | SA106C55 | 5*5ml |   |
| PreCap CHT I 60               | SA111C11 | 1*1ml |   |
|                               | SA111C15 | 1*5ml |   |
|                               | SA111C51 | 5*1ml |   |
|                               | SA111C55 | 5*5ml |   |

## Magnetic Beads


Magarose beads for protein purification

Magpoly beads for immunoprecipitation

# 04

## 4.1 Magnetic Beads

Magarose beads are superparamagnetic with rapid magnetic response, abundant hydroxyl functional groups, and uniform particle sizes, making them valuable tools in medical and molecular biology research. With particle sizes ranging from 30 to 100 μm, they are suitable for the detection and purification of target proteins.



**Magnetic beads:** Ni NTA Magarose Beads  
**Sample:** His-tagged protein expressed by E. coli  
**Binding / Wash Buffer:** 50 mM NaH<sub>2</sub>PO<sub>4</sub>, 0.3M NaCl, 20 mM imidazole, pH 8.0  
**Elution Buffer:** 50 mM NaH<sub>2</sub>PO<sub>4</sub>, 0.3 mM NaCl, 250 mM imidazole, pH 8.0

Fig.4-1 SDS-PAGE analysis of the purification of His-tagged protein using Ni-NTA Magarose Beads.

### Product Information

| Products                              | Cat.No.  | Size  | Characteristics (/ml medium) |                               |                          |
|---------------------------------------|----------|-------|------------------------------|-------------------------------|--------------------------|
| Ni IDA Magarose Beads                 | SM001001 | 1ml   | >40 mg His-tagged protein    |                               |                          |
|                                       | SM001005 | 5ml   |                              |                               |                          |
|                                       | SM001025 | 25ml  |                              |                               |                          |
|                                       | SM001100 | 100ml |                              |                               |                          |
| Ni NTA Magarose Beads                 | SM00101L | 1L    |                              |                               |                          |
|                                       | SM008001 | 1ml   |                              |                               |                          |
|                                       | SM008005 | 5ml   |                              |                               |                          |
|                                       | SM008025 | 25ml  |                              |                               |                          |
| Ni Smart Magarose Beads               | SM008100 | 100ml |                              |                               |                          |
|                                       | SM00801L | 1L    |                              |                               |                          |
|                                       | SM025001 | 1ml   |                              |                               |                          |
|                                       | SM025005 | 5ml   |                              |                               |                          |
| Glutathione Magarose Beads            | SM025025 | 25ml  | 10 mg His-tagged protein     |                               |                          |
|                                       | SM025100 | 100ml |                              |                               |                          |
|                                       | SM02501L | 1L    |                              |                               |                          |
|                                       | SM002001 | 1ml   |                              |                               |                          |
| rProtein A Magarose Beads             | SM002005 | 5ml   |                              | 5-10 mg GST-tagged protein    |                          |
|                                       | SM002025 | 25ml  |                              |                               |                          |
|                                       | SM002100 | 100ml |                              |                               |                          |
|                                       | SM00201L | 1L    |                              |                               |                          |
| rProtein G Magarose Beads             | SM003001 | 1ml   |                              |                               | >40mg Rabbit IgG         |
|                                       | SM003005 | 5ml   |                              |                               |                          |
|                                       | SM003025 | 25ml  |                              |                               |                          |
|                                       | SM003100 | 100ml |                              |                               |                          |
| Anti-DYKDDDDK Magarose Beads          | SM00301L | 1L    | >30mg Goat IgG               |                               |                          |
|                                       | SM004001 | 1ml   |                              |                               |                          |
|                                       | SM004005 | 5ml   |                              |                               |                          |
|                                       | SM004025 | 25ml  |                              |                               |                          |
| Anti-GFP Magarose Beads               | SM004100 | 100ml |                              | >1 mg DYKDDDDK-tagged protein |                          |
|                                       | SM00401L | 1L    |                              |                               |                          |
|                                       | SM009001 | 1ml   |                              |                               |                          |
|                                       | SM009005 | 5ml   |                              |                               |                          |
| NHS-Activated Magarose Beads          | SM009025 | 25ml  |                              |                               | >1 mg GFP-tagged protein |
|                                       | SM009100 | 100ml |                              |                               |                          |
|                                       | SM00901L | 1L    |                              |                               |                          |
|                                       | SM038001 | 1ml   |                              |                               |                          |
| Oligo(dT) <sub>25</sub> Magpoly Beads | SM038005 | 5ml   | >10mg Rabbit IgG             |                               |                          |
|                                       | SM038025 | 25ml  |                              |                               |                          |
|                                       | SM038100 | 100ml |                              |                               |                          |
|                                       | SM03801L | 1L    |                              |                               |                          |
| Epoxy-Activated MagPoly Beads         | SM031001 | 1ml   |                              |                               |                          |
|                                       | SM031005 | 5ml   |                              |                               |                          |
|                                       | SM031025 | 25ml  |                              |                               |                          |
|                                       | SM003100 | 100ml |                              |                               |                          |
|                                       | SM03101L | 1L    |                              |                               |                          |

## 4.2 Magpoly Beads for Immunoprecipitation

Magpoly beads designed for immunoprecipitation are coupled with a high density of ligands, making them have outstanding antibody binding capacity and low nonspecific adsorption. These beads can be used for one-step antibody purification from sample mixtures. Magpoly beads are widely used in immunoprecipitation, co-immunoprecipitation, and immunoadsorption applications.





Fig. 4-2 Magpoly Beads for Immunoprecipitation

### Product Information

| Products                              | Cat.No.  | Size  | Characteristics  |
|---------------------------------------|----------|-------|--|
| rProtein A MagPoly Beads              | SM037001 | 1ml   | >50 μg hlgG/mg Beads   |
|                                       | SM037005 | 5ml   |  |
|                                       | SM037010 | 10ml  |  |
|                                       | SM037050 | 50ml  |  |
|                                       | SM037100 | 100ml |  |
|                                       | SM037500 | 500ml |  |
| rProtein A/G MagPoly Beads            | SM03701L | 1L    | >50μg Rabbit IgG/mg Beads  |
|                                       | SM015001 | 1ml   |  |
|                                       | SM015005 | 5ml   |  |
|                                       | SM015010 | 10ml  |  |
|                                       | SM015050 | 50ml  |  |
|                                       | SM015100 | 100ml |  |
| Streptavidin MagPoly Beads            | SM015500 | 500ml | Intended for immobilization of biotin and capture of biotinylated substances including proteins, antibodies, and related biomolecules. |
|                                       | SM017001 | 1ml   |  |
|                                       | SM017005 | 5ml   |  |
|                                       | SM017010 | 10ml  |  |
|                                       | SM017050 | 50ml  |  |
|                                       | SM017100 | 100ml |  |
| NHS-Activated MagPoly Beads           | SM017500 | 500ml | It is pre-activated magpoly beads and can be directly used for the coupling of amino containing proteins or peptides.                  |
|                                       | SM01701L | 1L    |  |
|                                       | SM030001 | 1ml   |  |
|                                       | SM030005 | 5ml   |  |
|                                       | SM030010 | 10ml  |  |
|                                       | SM030050 | 50ml  |  |
| Oligo(dT) <sub>25</sub> Magpoly Beads | SM030100 | 100ml | Designed for mRNA isolation.   |
|                                       | SM030500 | 500ml |  |
|                                       | SM03001L | 1L    |  |
|                                       | SM080001 | 1ml   |  |
|                                       | SM080005 | 5ml   |  |
|                                       | SM080010 | 10ml  |  |
| Epoxy-Activated MagPoly Beads         | SM080050 | 50ml  | It is the superparamagnetic polymer beads with epoxy groups on the surface.  |
|                                       | SM080100 | 100ml |  |
|                                       | SM080500 | 500ml |  |
|                                       | SM08001L | 1L    |  |
|                                       | SM020001 | 1ml   |  |
|                                       | SM020005 | 5ml   |  |

## Instruments and Consumables

### High-throughput purification equipment

#### Protein chromatography system

#### Gravity Empty Column

#### Single-layer glass column


#### Double-layer column (SXX series)

# 05

## 5.1 High-throughput purification equipment



**MP10 (ten-channel protein purification system)** is suitable for process development and preparation at laboratory and industrial scales, ranging from analytical and small-scale preparation to pilot-scale operations. It utilizes chromatographic techniques such as gel filtration, ion exchange, affinity chromatography, and hydrophobic interaction chromatography to enable the separation, purification, and preparation of biomolecules including proteins, antibodies, peptides, polysaccharides, nucleic acids, as well as bioactive components from herbal medicines and natural products.



**Purifier 24** is designed for the automated transfer and processing of magnetic particles in a 24-well microplate format. Each well can accommodate a maximum volume of 4 mL, enabling optimal purification of free DNA, antibodies, and antigens. This system offers enhanced flexibility for fast and reproducible results across a wide range of applications.



**Purifier 32** uses magnetic bead separation technology. Unlike other methods that involve moving liquids, this system transfers magnetic beads between deep-well plates containing specific reagents. Magnetic rods equipped with disposable sleeves facilitate the collection, release, transfer, and incubation of magnetic beads by separating the rods from the sleeves. It can extract high-quality nucleic acids from samples such as saliva, swabs, dried blood spots, whole blood, serum/plasma, animal and plant tissues.

### Product Information

| abbreviation | Cat.No. |
|--------------|---------|
| MP 10        | E011    |
| Purifier 96  | P961001 |
| Purifier 32  | P321001 |
| Purifier 24  | P241001 |

## 5.2 Protein chromatography system



**Smart GO™** is a compact liquid chromatography system with highly integrated piping design and modular layout, minimizing system dead volume and optimizing experimental space. The Vaticy operating software, with built-in parameter data of various Smart-Lifesciences pre-packed columns, making it faster and easier for users to start and operate. It is mainly used for laboratory-scale process research, and can be used to purify proteins, antibodies, vaccines, peptides, viruses, and other macromolecules, ranging from micrograms to grams in size.



**Smart PUR™** is a precise and stable chromatography system with flexible expansion configurations, supporting the purification of recombinant proteins, antibodies, vaccines, nucleic acids, diagnostic materials, and other target products, from microgram to gram quantities, making it suitable for the diverse experimental studies. The Vaticy operating software, with built-in parameter data of various Smart-Lifesciences pre-packed columns, making it faster and easier for users to start and operate, and is the best partner for process development.




**SmartPurify Lite™** is an efficient and intelligent protein purification device designed for applications in biopharmaceuticals, research institutions, and related fields. Featuring a modular and integrated architecture, the system optimizes laboratory space utilization. It is equipped with a 12-inch touch screen and intelligent guidance functions, supporting custom programming and one-click operation to streamline chromatography workflows. With high-precision dual-plunger pumps, a UV detector, and a real-time pressure monitoring system, it enables precise control of separation conditions, providing reliable experimental support for applications such as recombinant protein expression and antibody purification.

### Product Information

| Products         | Cat.No. | Size |
|------------------|---------|------|
| Smart GO 150     | E015    | unit |
| Smart GO 100     | E016    | unit |
| Smart GO 35      | E017    | unit |
| Smart PUR 150    | E012    | unit |
| Smart PUR 100    | E013    | unit |
| Smart PUR 35     | E014    | unit |
| SmartPurify Lite | E021    | unit |

## 5.3 Gravity Empty Column

| Empty column of gravity chromatogra- |  |        |        |         |         |         |         |
|--------------------------------------|--|--------|--------|---------|---------|---------|---------|
| Cat.No.                              | SLM001                                     | SLM004 | SLM005 | SLM006  | SLM007  | SLM025  | SLM009  |
| Size (ml)                            | 12 ml                                      | 1ml    | 3 ml   | 6 ml    | 30 ml   | 60 ml   | 300 ml  |
| Inner diameter (mm)                  | 15.6 mm                                    | 5.6 mm | 8.9 mm | 12.7 mm | 22.5 mm | 26.2 mm | 49.4 mm |
| height (mm)                          | 83 mm                                      | 55 mm  | 63 mm  | 66 mm   | 101 mm  | 134 mm  | 180 mm  |
| volume (ml)                          | 12 ml                                      | 1 ml   | 3 ml   | 6 ml    | 30 ml   | 60 ml   | 300 ml  |
| Material (column tube)               | polypropylene                              |        |        |         |         |         |         |
| Material (filter membrane)           | polyethylene                               |        |        |         |         |         |         |
| Pore size of filter                  | 50µm                                       |        |        |         |         |         |         |
| Chemical stability                   | It is stable in all commonly used reagents |        |        |         |         |         |         |



## 5.4 Single-layer glass column

**Column Structure:** The empty column, upper flange, adapter (axial expansion), sealing ring, lower flange, inlet and outlet.

**Material:** High borosilicate glass.

**Maximum Pressure:** 3-5 bar.

**Chemical Stability:** Stable in the presence of acids, alkalis, and organic solvents (Should not be operated with solvents such as acetone, methane, ammonium hydroxide, or dichloromethane for extended periods).

**Axial Expansion:** The axial expansion feature allows for adjustments in packing height, thereby enhancing separation efficiency.

**Operating Temperature:** 0 °C – 60 °C.

**Filter:** 50 mesh –500 mesh, depending on the size of the medium, effectively blocking general dextran gel, agarose gel, and macroporous resin.

**Cleaning:** Soapy water and detergent are suitable for cleaning the columns. If the column is contaminated by protein, it can be cleaned with detergent.


**Sterilization:** Under 120 °C, autoclave for 20 minutes.

### Product Information

| Single-layer Column (S Series) |         |                    |                |         |                    |
|--------------------------------|---------|--------------------|----------------|---------|--------------------|
| inner diameter                 | Cat.No. | Size (I.D. × L mm) | inner diameter | Cat.No. | Size (I.D. × L mm) |
| 10mm                           | S1020   | 10*200             | 15mm           | S1520   | 15*200             |
|                                | S1030   | 10*300             |                | S1530   | 15*300             |
|                                | S1040   | 10*400             |                | S1540   | 15*400             |
|                                | S1050   | 10*500             |                | S1550   | 15*500             |
|                                | S1060   | 10*600             |                | S1560   | 15*600             |
|                                | S1075   | 10*750             |                | S1575   | 15*750             |
|                                | S10100  | 10*1000            |                | S15100  | 15*1000            |
| 25mm                           | S2520   | 25*200             | 35mm           | S3520   | 50*200             |
|                                | S2530   | 25*300             |                | S3530   | 50*300             |
|                                | S2540   | 25*400             |                | S3540   | 50*400             |
|                                | S2550   | 25*500             |                | S3550   | 50*500             |
|                                | S2560   | 25*600             |                | S3560   | 50*600             |
|                                | S2575   | 25*750             |                | S3575   | 50*750             |
|                                | S25100  | 25*1000            |                | S35100  | 50*1000            |
| 50mm                           | S5020   | 50*200             | 75mm           | S7530   | 75*300             |
|                                | S5030   | 50*300             |                | S7540   | 75*400             |
|                                | S5040   | 50*400             |                | S7550   | 75*500             |
|                                | S5050   | 50*500             |                | S7560   | 75*600             |
|                                | S5060   | 50*600             |                | S7575   | 75*750             |
|                                | S5075   | 50*750             |                | S75100  | 75*1000            |
|                                | S50100  | 50*1000            |                |         |                    |

Product Information

| Single-layer Column (S Series) |         |                    |
|--------------------------------|---------|--------------------|
| inner diameter                 | Cat.No. | Size (I.D. × L mm) |
| 100mm                          | S10030  | 100*300            |
|                                | S10040  | 100*400            |
|                                | S10050  | 100*500            |
|                                | S10060  | 100*600            |
|                                | S10075  | 100*750            |
|                                | S100100 | 100*1000           |



## 5.5 Double-layer Column (SXX Series)

### 01 Uniform Flow Rate

Plunger design for ensuring homogeneous buffer distribution and reproducible packing.

### 02 Temperature Stability

Adapter with QuickLock mechanism for enhanced column handling and cleaning.

### 03 Good Column Efficiency

Fittings for direct connection to chromatography systems.

### 04 High Chemical Tolerance

Cooling jacket for sensitive purifications.

Product Information

| Double-layer Column (SXX Series) |                    |                |                     |                     |
|----------------------------------|--------------------|----------------|---------------------|---------------------|
| Cat.No.                          | Size (I.D. × L mm) | Pressure (bar) | Filling volume (ml) | Loading height (mm) |
| SXX 16/20                        | 16*200             | 5              | 4-33                | 20-165              |
| SXX 16/30                        | 16*300             | 5              | 24-53               | 120-265             |
| SXX 16/40                        | 16*400             | 5              | 44-73               | 220-365             |
| SXX 16/50                        | 16*500             | 5              | 64-93               | 320-465             |
| SXX 16/60                        | 16*600             | 5              | 84-113              | 420-565             |
| SXX 16/70                        | 16*700             | 5              | 104-133             | 520-665             |
| SXX 16/80                        | 16*800             | 5              | 124-153             | 620-765             |
| SXX 16/90                        | 16*900             | 5              | 144-173             | 720-865             |
| SXX 16/100                       | 16*1000            | 5              | 164-193             | 820-965             |
| SXX 26/20                        | 26*200             | 5              | 10.5-86             | 20-165              |
| SXX 26/30                        | 26*300             | 5              | 63-140              | 120-265             |
| SXX 26/40                        | 26*400             | 5              | 115-190             | 220-365             |
| SXX 26/50                        | 26*500             | 5              | 167-242             | 320-465             |
| SXX 26/60                        | 26*600             | 5              | 220-295             | 420-565             |
| SXX 26/70                        | 26*700             | 5              | 273-350             | 520-665             |
| SXX 26/80                        | 26*800             | 5              | 325-400             | 620-765             |
| SXX 26/90                        | 26*900             | 5              | 377-455             | 720-865             |
| SXX 26/100                       | 26*1000            | 5              | 430-508             | 820-965             |
| SXX 50/20                        | 50*200             | 3              | 30-330              | 15-165              |
| SXX 50/30                        | 50*300             | 3              | 230-530             | 115-265             |
| SXX 50/40                        | 50*400             | 3              | 430-730             | 215-365             |
| SXX 50/50                        | 50*500             | 3              | 630-930             | 315-465             |
| SXX 50/60                        | 50*600             | 3              | 830-1130            | 415-565             |
| SXX 50/70                        | 50*700             | 3              | 1030-1330           | 515-665             |
| SXX 50/80                        | 50*800             | 3              | 1230-1530           | 615-765             |
| SXX 50/90                        | 50*900             | 3              | 1430-1730           | 715-865             |
| SXX 50/100                       | 50*1000            | 3              | 1630-1930           | 815-965             |

| Double-layer Column (SXX Series) |                    |                |                     |                     |
|----------------------------------|--------------------|----------------|---------------------|---------------------|
| Cat.No.                          | Size (I.D. × L mm) | Pressure (bar) | Filling volume (ml) | Loading height (mm) |
| SXX16/20D                        | 16*200             | 5              | 0-33                | 20-165              |
| SXX16/30D                        | 16*300             | 5              | 0-53                | 0-265               |
| SXX16/40D                        | 16*400             | 5              | 15-73               | 75-365              |
| SXX16/50D                        | 16*500             | 5              | 35-93               | 175-465             |
| SXX16/60D                        | 16*600             | 5              | 55-113              | 275-565             |
| SXX16/70D                        | 16*700             | 5              | 75-133              | 375-665             |
| SXX16/80D                        | 16*800             | 5              | 95-153              | 475-765             |
| SXX16/90D                        | 16*900             | 5              | 115-173             | 575-865             |
| SXX16/100D                       | 16*1000            | 5              | 135-193             | 675-965             |
| SXX26/20D                        | 26*200             | 5              | 0-86                | 0-165               |
| SXX26/30D                        | 26*300             | 5              | 0-140               | 0-265               |
| SXX26/40D                        | 26*400             | 5              | 38-190              | 75-365              |
| SXX26/50D                        | 26*500             | 5              | 87-242              | 175-465             |
| SXX26/60D                        | 26*600             | 5              | 138-295             | 275-565             |
| SXX26/70D                        | 26*700             | 5              | 190-350             | 375-665             |
| SXX26/80D                        | 26*800             | 5              | 240-400             | 475-765             |
| SXX26/90D                        | 26*900             | 5              | 290-455             | 575-865             |
| SXX26/100D                       | 26*1000            | 5              | 340-508             | 675-965             |
| SXX50/20D                        | 50*200             | 3              | 0-330               | 0-165               |
| SXX50/30D                        | 50*300             | 3              | 0-530               | 0-265               |
| SXX50/40D                        | 50*400             | 3              | 160-730             | 80-365              |
| SXX50/50D                        | 50*500             | 3              | 360-930             | 180-465             |
| SXX50/60D                        | 50*600             | 3              | 560-1130            | 280-565             |
| SXX50/70D                        | 50*700             | 3              | 760-1330            | 380-665             |
| SXX50/80D                        | 50*800             | 3              | 960-1530            | 480-765             |
| SXX50/90D                        | 50*900             | 3              | 1160-1730           | 580-865             |
| SXX50/100D                       | 50*1000            | 3              | 1360-1930           | 680-965             |

Product Information

| SXX Series Accessories |         |                                      |
|------------------------|---------|--------------------------------------|
| Products               | Cat.No. | Notes                                |
| Packing reservoir      | SXX1601 | Compatible with SXX16 series columns |
| adapter                | SXX1602 |                                      |
| bottom plunger         | SXX1603 |                                      |
| filter                 | SXX1604 |                                      |
| filter-holder          | SXX1605 |                                      |
| O-ring                 | SXX1606 |                                      |
| Packing connector      | SXX1607 | Compatible with SXX26 series columns |
| Packing reservoir      | SXX2601 |                                      |
| adapter                | SXX2602 |                                      |
| bottom plunger         | SXX2603 |                                      |
| filter                 | SXX2604 |                                      |
| filter-holder          | SXX2605 |                                      |
| O-ring                 | SXX2606 |                                      |
| Packing connector      | SXX2607 |                                      |

| Products          | Cat.No. | Notes                                |
|-------------------|---------|--------------------------------------|
| Packing reservoir | SXK5001 | Compatible with SXK50 series columns |
| adapter           | SXK5002 |                                      |
| bottom plunger    | SXK5003 |                                      |
| filter            | SXK5004 |                                      |
| filter-holder     | SXK5005 |                                      |
| O-ring            | SXK5006 |                                      |
| Packing connector | SXK5007 |                                      |

